



EDINBURGH
INSTRUMENTS

FLS1000

User Guide



Copyrights

Copyright © 2020 Edinburgh Instruments Ltd. All rights reserved.

No part of this publication may be reproduced, stored in a retrieval system, published on the internet or transmitted in any form or any means electronic or mechanical, including photocopying and recording for any purpose other than the purchaser's personal use without written permission.

Information in this document is subject to change without notice. While every effort has been taken to make sure that the information given is correct, Edinburgh Instruments Limited cannot accept any liability for errors and omissions or their consequences. The software described in this document is furnished under a license agreement or nondisclosure agreement. The software may be used or copied only in accordance with the terms of those agreements.

Edinburgh Instruments Ltd.

2 Bain Square

Kirkton Campus

Livingston

EH54 7DQ

UK

www.edinst.com

FLS1000 User Guide

Issue No. 3

Issue Date: 30 January 2020

© Edinburgh Instruments Ltd
2 Bain Square, Kirkton Campus, Livingston, EH54 7DQ UK
Phone: +44 (0) 1506 425300 Fax: +44 (0) 1506 425320
Email: eai.service@edinst.com
<http://www.edinst.com>

Contents

OVERVIEW	1
1. INTRODUCTION TO THE FLS1000 SERIES OF FLUORESCENCE SPECTROMETERS	3
1.1. FLS1000 SERIES MODELS	4
1.2. SPECTROMETER UPGRADES	11
2. SPECTROMETER COMPONENTS	13
2.1. LIGHT SOURCES	13
2.1.1. <i>Xe2 Continuous xenon Lamp.....</i>	<i>13</i>
2.1.2. <i>µF2 Xenon Flashlamp.....</i>	<i>15</i>
2.1.3. <i>Picosecond Pulsed Diode Laser (EPL) and Picosecond Pulsed LED (EPLD).....</i>	<i>17</i>
2.1.4. <i>Other Light Sources.....</i>	<i>17</i>
2.2. MONOCHROMATORS.....	18
2.2.1. <i>Excitation Monochromator.....</i>	<i>18</i>
2.2.2. <i>Emission Monochromator</i>	<i>20</i>
2.2.3. <i>Monochromator Gratings.....</i>	<i>22</i>
2.2.4. <i>Slits</i>	<i>22</i>
2.2.5. <i>Shutters.....</i>	<i>22</i>
2.2.6. <i>Filter Wheel.....</i>	<i>23</i>
2.3. SAMPLE CHAMBER	24
2.3.1. <i>Standard Lens Optics</i>	<i>26</i>
2.3.2. <i>Mirror Optics and Other Beam Steering Options.....</i>	<i>26</i>
2.4. SAMPLE HOLDERS	27
2.4.1. <i>Standard Single Cuvette Holder</i>	<i>27</i>
2.4.2. <i>Other Sample Holders and Sample Positioning Options</i>	<i>27</i>
2.5. DETECTORS	28
2.5.1. <i>Red Sensitive Photomultiplier in Cooled Housing</i>	<i>28</i>
2.5.2. <i>Other Detectors</i>	<i>29</i>
3. INTRODUCTION TO FLUORACLE® SOFTWARE.....	30
3.1. STARTING FLUORACLE	30
3.2. GENERAL FEATURES OF THE FLUORACLE MAIN WINDOW	30
3.2.1. <i>Menu Bar</i>	<i>30</i>
3.2.2. <i>Tool Bar.....</i>	<i>30</i>
3.2.3. <i>Scan Display Area</i>	<i>32</i>
3.2.4. <i>Status Bar</i>	<i>32</i>
3.2.5. <i>Signal Rate Window.....</i>	<i>32</i>
3.3. SCAN TYPES AND SCAN PROPERTIES	34
3.3.1. <i>Spectral Data</i>	<i>34</i>
3.3.2. <i>Time-Resolved Data.....</i>	<i>36</i>
3.4. GRAPHICAL PRESENTATIONS AND PLOT OPTIONS	37
3.4.1. <i>2D Data Presentation</i>	<i>37</i>
3.4.2. <i>3D Data Presentation</i>	<i>40</i>
3.4.3. <i>Contour Plots/Color Maps</i>	<i>41</i>
3.4.4. <i>View Numerical Data.....</i>	<i>43</i>
3.5. DATA OPERATION WITH FLUORACLE.....	44
3.6. ANALYSIS WITH FLUORACLE.....	45
3.7. SOFTWARE AND HARDWARE OPTIONS.....	46
4. STARTING UP, SETTING MEASUREMENT AND ANALYSIS PREFERENCES AND SHUTTING DOWN	47
4.1. START-UP PROCEDURE	47
4.2. SETTING UP OPTIONS AND PREFERENCES	48
4.2.1. <i>Allocate Correction File(s) to the Light Path(s)</i>	<i>48</i>

4.2.2.	Select Suitable Window Options	48
4.2.3.	Selecting Measurement Options	49
4.2.4.	Selecting signals to be displayed on Signal Rate dialogue	49
4.2.5.	Enabling Special Analysis Options	50
4.2.6.	Setting Options for Data Export	50
4.2.7.	Selecting the Sample Holder	51
4.3.	SHUTTING DOWN THE SPECTROMETER	52
5.	MAKING SPECTRAL MEASUREMENTS	53
5.1.	PREPARING FOR A MEASUREMENT – THE SIGNAL RATE DIALOGUE	54
5.2.	EMISSION SCANS	56
5.3.	EXCITATION SCANS	57
5.4.	SYNCHRONOUS SCANS	58
5.6.	EMISSION MAPS	60
5.7.	SYNCHRONOUS MAPS	62
5.8.	MULTIPLE SPECTRAL MEASUREMENTS	63
	TIPS FOR MAKING GOOD SPECTRAL MEASUREMENTS	65
6.	MAKING TIME-RESOLVED MEASUREMENTS	69
6.1.	KINETIC MEASUREMENTS	72
6.2.	MANUAL TIME-RESOLVED MEASUREMENTS – MCS	74
6.2.1.	<i>µF Lamp Setup</i>	74
6.2.2.	<i>Measurement of a Phosphorescence Decay</i>	74
6.2.3.	<i>Measurement of the Instrument Response Function</i>	77
6.3.	MANUAL TIME-RESOLVED MEASUREMENTS – TCSPC	78
6.3.1.	<i>EPL Setup</i>	78
6.3.2.	<i>Measurement of a Fluorescence Decay</i>	78
6.3.3.	<i>Measurement of the Instrument Response Function</i>	81
6.4.	MULTIPLE TIME-RESOLVED MEASUREMENTS	82
6.5.	TRES-MAP MEASUREMENTS	83
6.5.1.	<i>Excitation TRES-Map Measurements</i>	83
6.5.2.	<i>Emission TRES-Map Measurements</i>	83
6.6.	TIPS FOR MAKING GOOD LIFETIME MEASUREMENTS	84
7.	RE-RUN AND BATCH MEASUREMENTS	88
7.1.	RE-RUN A SPECTRAL MEASUREMENT	88
7.2.	BATCH MEASUREMENTS	88
8.	DATA OPERATIONS AND ANALYSIS ON SPECTRAL DATA	90
8.1.	ADD, SUBTRACT, MULTIPLY, DIVIDE	90
8.2.	APPEND	91
8.3.	SCALE, NORMALISE, SUBTRACT BASELINE, CROP RANGE	92
8.4.	SMOOTH	94
8.5.	DIFFERENTIATE AND INTEGRATE	95
8.6.	SPECTRAL CORRECTION	96
8.7.	CHROMATICITY PLOTS	97
9.	DATA OPERATIONS AND ANALYSIS ON TIME-RESOLVED DATA	99
9.1.	ADD, SUBTRACT, MULTIPLY, DIVIDE	99
9.2.	SCALE, NORMALISE, SUBTRACT BASELINE, CROP RANGE	99
9.3.	ANALYSIS OF LUMINESCENCE KINETICS	99
9.3.1.	<i>Tail Fit Analysis</i>	100
9.3.2.	<i>Reconvolution Fit Analysis</i>	104
9.4.	FIT QUALITY PARAMETERS	108
9.4.1.	<i>The Reduced Chi-Square</i>	108
9.4.2.	<i>The Residual Data</i>	108

9.4.3. <i>The Autocorrelation Data</i>	109
9.4.4. <i>The Durbin-Watson Parameter</i>	109
9.5. TRES DATA SLICING	110
FORMULAS AND DEFINITIONS	113
GLOSSARY	114
INDEX	115
TABLE OF FIGURES	117
TABLE OF PROCEDURES	121

Overview

The FLS1000 Series of Fluorescence Spectrometers enable steady state, fluorescence lifetime and phosphorescence lifetime measurements.

Read this document if you plan to use the spectrometer to make measurements or use the software to analyse your measurement results. It describes the spectrometer components, the software you use to control the spectrometer and how to make spectral and lifetime measurements.

This User Guide covers the operation of the standard components. Making measurements are explained in chapters 5 and 6. The data operation and analysis procedures are explained in chapters 8 and 9.

More Information

For details on how to install an FLS1000 spectrometer, for service procedures and for validation procedures refer to the *FLS1000 Series Installation and Service Guide*.

Special Accessories and Hardware options will be supplied with their own reference guides. For details on these items refer to the specific reference guide.

1. Introduction to the FLS1000 Series of Fluorescence Spectrometers

The FLS1000 Series is a series of modular software-controlled spectrometers for the acquisition of steady state and time-resolved photoluminescence. They can measure:

- Fluorescence and phosphorescence spectra
- Fluorescence and phosphorescence lifetimes

Data can be obtained over a wide spectral range, from the ultraviolet to the near-infrared, with single photon counting sensitivity. The spectrometer series combines superior sensitivity with high spectral and temporal resolution.

The FLS1000 Series of Fluorescence Spectrometers is ideal for demanding applications in such areas as:

- Photophysics
- Photochemistry
- Biophysics
- Materials research

An FLS1000 Series spectrometer has the following main components:

- *Light Source* – an excitation source to generate photoluminescence from the sample
- *Excitation Monochromator* – for the selection of specific monochromatic light from the full spectrum of the excitation source
- *Sample Chamber* – for optical components and sample holders
- *Emission Monochromator* – for the selection of specific monochromatic light from the full spectrum of the sample emission
- *Detector* – detects the photon flux of the sample emission at the selected wavelength
- *Dedicated PC* – is pre-installed with spectrometer operating software Fluoracle.

Depending on the type of measurements you want to make, different light sources are used together with different data acquisition techniques:

Measurement Type	Standard Light Source	Data Acquisition Technique
Spectral measurements	Continuous xenon lamp (Xe2)	Single Photon Counting
Time-resolved measurements in the microsecond to second range	Microsecond flashlamp (μ F2)	Multi-channel scaling (MCS) for time-resolved photon counting
Time-resolved measurements in the pico- to nano-second range	Picosecond pulsed lasers EPL or EPLED, nF920 nanosecond flashlamp as an option	Time-Correlated Single Photon Counting (TCSPC)

The following section (chapter 1.1) introduces each of the ten models in the FLS1000 Series of fluorescence spectrometers. Your own spectrometer may differ from those shown in respect to light sources, detectors and optional accessories. Upgrade options are listed in chapter 1.2 on page 11.

1.1. FLS1000 Series Models

The ten FLS1000 Series models are listed below.

FLS1000-SS-s and FLS1000-DD-s

The FLS1000-SS-s and FLS1000-DD-s are spectrometers for the acquisition of fluorescence spectra and fluorescence kinetic measurements. The FLS1000-DD-s is the world's most sensitive steady-state fluorescence spectrometer, with a specified signal-to-noise ratio of the Water Raman signal of 30,000:1. The signal-to-noise ratio of the FLS1000-SS-s specifies at 20,000:1, while there are also FLS1000-SD-s and FLS1000-DS-s models with signal-to-noise ratio of 25,000:1.

Figure 1-1a shows the layout of the standard version of the FLS1000-SS-s with single grating monochromators; Figure 1-1b shows the FLS1000-DD-s with double grating monochromators.

Component	Standard
Light Source	Xe2 continuous xenon lamp
Detector	Single Photon Counting PMT
Data Acquisition Module	CB1
Acquisition Technique	Single Photon Counting (spectral scanning, kinetic measurements)




Figure 1-1: General layout of, a) the FLS1000-SS-s and b) the FLS1000-DD-s

FLS1000-xS-t and FLS1000-SS-t

The FLS1000-xS-t and FLS1000-SS-t are dedicated spectrometer for the acquisition of fluorescence decay kinetics the time range from picoseconds to microseconds, based on the technique of Time-Correlated Single Photon Counting (TCSPC).

Component	Standard
Light Source	EPLs and EPLEDs (nF920 nanosecond flashlamp as an option)
Detector	Single Photon Counting PMT
Data Acquisition Card	TCC2
Acquisition Technique	TCSPC (fluorescence decay acquisitions)




Figure 1-2: General layout of a) the FLS1000-xS-t with EPL/EPLED, b) the FLS1000-SS-t with nF920.

FLS1000-SS-m

The FLS1000-SS-m is a dedicated spectrometer for the acquisition of phosphorescence decay kinetics in the time range from microseconds to seconds, based on the technique of multi-channel scaling (MCS).

Component	Standard
Light Source	μ F2 microsecond flashlamp
Detector	Single Photon Counting PMT
Data Acquisition Card	CB1
Acquisition Technique	MCS (phosphorescence decay acquisitions)




Figure 1-3: General layout of an FLS1000-SS-m

FLS1000-SS-sm

The FLS1000-SS-sm is a steady state and phosphorescence lifetime spectrometer. It can be used to measure fluorescence and phosphorescence spectra, as well as phosphorescence decay kinetics.

Component	Standard
Light Sources	Xe2 continuous xenon lamp, μ F1 microsecond flashlamp
Detector	Single Photon Counting PMT
Data Acquisition Card	CB1
Acquisition Technique	MCS (spectral scanning and phosphorescence decay acquisitions)




Figure 1-4: General layout of an FLS1000-SS-sm

FLS1000-SS-tm

The FLS1000-SS-tm is a dedicated instrument for the acquisition of fluorescence and phosphorescence kinetics throughout all time ranges (picoseconds to seconds).

Component	Standard
Light Sources	μ F2 microsecond flashlamp, EPLs/EPLEDs (nF920 nanosecond flashlamp as an option)
Detector	Single Photon Counting PMT
Data Acquisition Card	CB1, TCC2
Acquisition Technique	MCS and TCSPC (fluorescence and phosphorescence decay acquisitions)




Figure 1-5: General layout of an FLS1000-SS-tm

FLS1000-SS-st

The FLS1000-SS-st is a combined steady state and time-resolved fluorescence spectrometer.

Component	Standard
Light Sources	Xe2 continuous xenon lamp, EPL/EPLEDs (nF920 nanosecond flashlamp as an option)
Detector	Single Photon Counting PMT
Data Acquisition Card	CB1, TCC2
Acquisition Technique	Single Photon Counting, MCS and TCSPC (spectral scanning, kinetic measurements, fluorescence decay acquisitions)




Figure 1-6: General layout of an FLS1000-SS-st

FLS1000-SS-stm and FLS1000-DD-stm

The FLS1000-SS-stm and FLS1000-DD-stm are complete photoluminescence laboratories in a single instrument. They combine all the features of the abovementioned models as detailed below.

Figure 1-7a shows the layout of the standard version of the FLS1000-SS-stm with single grating monochromators, Figure 1-7b shows the FLS1000-DD-stm with the popular upgrade of double grating monochromators.

Component	Standard
Light Sources	Xe2 continuous xenon lamp, μ F2 microsecond flashlamp, EPLs and EPLEDs (nF920 nanosecond flashlamp as an option)
Detector	Single Photon Counting PMT
Data Acquisition Card	CB1, TCC2
Acquisition Technique	Single Photon Counting, MCS and TCSPC (spectral scanning, kinetic measurements, fluorescence and phosphorescence decay acquisitions)




Figure 1-7: General layout of a) an FLS1000-SS-stm, b) an FLS1000-DD-stm.

1.2. Spectrometer Upgrades

The fluorescence spectrometers listed in Section 1.1 are standard configurations. Various upgrades, modifications, and enhancements are also available, giving the FLS1000 the flexibility to be used with a wide range of specific applications and sample requirements.

Common upgrades and enhancements include:

- Double monochromators with increased stray light performance for better measurements of highly scattering samples
- Additional emission monochromators to form "T" or "X" geometries, enabling more than two detectors to be permanently installed (or more than three when using double monochromators)
- Polarisers for automatic acquisition of fluorescence anisotropy data
- Gated photomultipliers for better discrimination between (long) phosphorescence decay and strong (short) fluorescence and/or strong sample scattering
- Near Infrared detectors

Additionally, a range of sample holders and sample cooling options are available:

- Cuvette holders with options for water circulation (for cooling or heating) with a temperature probe
- Magnetic stirrers
- Front face sample holders on slides, rotational stages, X-Y stages
- Thermo-electrically temperature-controlled cuvette holder with magnetic stirrer option
- EPR dewar option for cooling at liquid nitrogen temperature
- Cryostat options (Helium or Nitrogen)
- Multi-position sample holders, with water circulation or thermo-electric based temperature control
- Plate Reader attachment
- Fibre attachments
- Fibre couplers to microscopes
- Integration Sphere accessory for the measurement of absolute quantum yield of fluorescent or phosphorescent samples
- Absorption accessory
- X-ray excitation

For more details on these options, please contact Edinburgh Instruments.

2. Spectrometer Components

The following sections describe the FLS1000 Series standard components. For details on other components, please contact Edinburgh Instruments.

2.1. Light Sources

There are three standard light sources in the FLS1000 spectrometer series:

- A 450 W continuous xenon arc lamp for continuous sample excitation.
- A low repetition rate pulsed 60W xenon flashlamp for pulsed sample excitation in multichannel scaling experiments.
- EPLs (picosecond pulsed diode lasers) or EPLEDs (picosecond pulsed LEDs) for Time-Correlated Single Photon Counting experiments.

2.1.1. Xe2 Continuous xenon Lamp

The standard light source for steady state applications is a 450 W xenon arc lamp, as shown in Figure 2-1. It produces "white" light composed of a continuum superimposed with narrow peaks. The optimal spectral range extends from below 250 nm to more than 1000 nm, shown in Figure 2-2.




Figure 2-1: Xe2 Xenon Lamp

The Xe2 is an integrated module, comprising focusing optics, power supply, ignitor and control electronics. The emission of the arc is focussed onto the slit of the excitation monochromator using an ellipsoidal mirror. The lamp has a removable service panel to access the bulb for alignment or replacement. The display panel at the front of the lamp indicates the operational status, count, voltage, power status and elapsed operation time.




Figure 2-2: Spectral Output of the 450 W Xenon Arc Lamp

The xenon lamp can be fitted with an ozone-free bulb or with an ozone generating bulb (optional). When fitted with an ozone generating bulb, the spectral output in the UV spectral range is enhanced and the lower wavelength limit shifts from about 250 nm to 230 nm. Further enhancement in the UV can be achieved by selecting special UV gratings for the excitation monochromator.

	<p>Caution Ozone</p>	<p>When working with ozone-generating bulbs, an air extraction system must be fitted to the lamp head.</p>
--	---------------------------------	---

Maintenance

The Xe2 is mostly maintenance-free. Optical alignment is made within the factory, optimal alignment will be checked during installation; see the *FLS1000 Series Installation and Service Guide* for details.

The xenon bulb has an average lifetime of 2000 hours (1000 hours guaranteed). For instructions on how to change and align the bulb, refer to the *FLS1000 Series Installation and Service Guide* for details.

2.1.2. μ F2 Xenon Flashlamp

The μ F2 is a compact 60 W xenon flashlamp, which is triggered by the spectrometer's controller. The output is focussed by an adjustable lens assembly onto one of the entrance slits in the spectrometer's excitation monochromator.



Figure 2-3: μ F2 Xenon Flashlamp

The μ F2 is optimised for multi-channel scaling (MCS) lifetime measurements. It has a narrow optical pulse width of $1\ \mu\text{s} - 2\ \mu\text{s}$ and a pulse repetition rate of between $0.1\ \text{Hz} - 100\ \text{Hz}$.

For some applications (for example, time-gated spectral scanning), the μ F2 lamp can be used as an excitation source for spectral scans. When operating in this mode your samples should have intrinsically long lifetimes since the spectrometer uses the highly sensitive single photon counting mode for spectral scanning. The number of instantaneous photons counted per dwell time is limited by a repetition rate of $100\ \text{Hz}$ and a typical pulse width of $1\ \mu\text{s} - 2\ \mu\text{s}$.

μF2 Output Spectrum

The typical output spectrum of the μF2 is shown in Figure 2-4. The microsecond pulsed xenon flashlamp has an optical output that is more structured with narrow peaks and has comparatively higher content of UV output when compared to the spectrum of a continuous xenon lamp (see Section 2.1.1 Xe2 Continuous xenon Lamp on page 13).




Figure 2-4: Spectral Output of the Xenon Flashlamp

μF2 Pulse Profile

Figure 2-5 shows the temporal output of the μF2 flashlamp, the typical pulse width (FWHM) is 2 μs. Note that the shape and the tail of the pulse is not equal for all wavelengths.




Figure 2-5: Temporal Output of μF2 Flashlamp at 300 nm (blue), 450 nm (green) and 600 nm (red).

2.1.3. Picosecond Pulsed Diode Laser (EPL) and Picosecond Pulsed LED (EPLD)

Pulsed diode lasers and LEDs are the most popular light sources for lifetime measurements in the nanosecond range for the FLS1000 Series of spectrometers. Edinburgh Instruments can supply a range of picoseconds pulsed laser diodes (EPLs) and LEDs (EPLDs) spanning the spectrum from 250 nm – 980 nm. These are connected to the spectrometer using a coupling flange dedicated to pulsed light sources. The coupling flange connects to the spare port of the sample chamber opposite the excitation arm. This flange also has an integrated adjustable neutral density filter.



Figure 2-6: Edinburgh Instruments Picosecond Pulsed Diode Laser

All picosecond pulsed light sources require a beam dump and a beam aperture for precise measurements in the lower picosecond region. The beam dump accepts the residual beam after passing through the sample and attenuates it. This significantly reduces back reflections from lens surfaces and other components in the sample chamber that can impair the time-resolved measurements.

In addition, a beam aperture can be inserted into the cuvette sample holder, limiting the field of view for the emission channel. This reduces the temporal spread that is caused by the propagation time of the laser beam travelling through the cuvette (~30 ps).

Beam dump and beam aperture are not usually required for EPLDs due to the lower power output compared to EPLs. However, they can significantly improve data quality when EPLs are used with fast detectors.

2.1.4. Other Light Sources

For details on other light sources such as NIR CW lasers Nd: YAG lasers, Supercontinua, please contact Edinburgh Instruments.

2.2. Monochromators

A *monochromator* transmits a selectable narrow band of wavelengths of light from a wider range of wavelengths available at the input. Most FLS1000 fluorescence spectrometers have at least two monochromators: an *excitation monochromator* and an *emission monochromator*. These can be either single-grating, double-grating or triple grating monochromators, depending on your configuration. The FLS1000 can be configured in a "T" geometry with additional emission monochromator. This enables several detectors to be permanently installed.

FLS1000 fluorescence spectrometers use symmetric Czerny-Turner type monochromators, specifically designed to meet the requirements of the modular concept of the spectrometer. The monochromators have two chambers:

- An upper optical chamber, coated with highly absorbing optical grade paint of less than 0.2% reflectance. It is fitted with baffles and shields to minimise reflections of scattered light on reflecting surfaces and is sealed against ambient light for single photon integrity.
- The lower chamber contains mechanical and electrical components

The focal length of the monochromators is 325 mm with an F-number of 4. They are fitted with one (single monochromators) or two (double monochromators) grating turrets, each one carrying up to three gratings. The monochromators also have several computer-controlled slits, a filter wheel, swing mirrors and shutters.

When you load the spectrometer software, the monochromators calibrate by finding defined positions for the grating turrets, the slit mechanism and the filter wheel.

Note: There are lower-priced versions of the monochromators available with reduced number of slits and without the filter wheel.

2.2.1. Excitation Monochromator

The single-grating excitation monochromator has two entrance ports for lamps and one exit port mounted to the spectrometer's sample chamber:




Figure 2-7: TGM325-X Single-Grating Excitation Monochromator

Entrance port B is fitted with a computer-controlled shutter. Although typically configured for the Xe2 continuous xenon lamp, it may also be used for the μ F2 microsecond flashlamp.




Figure 2-8: Schematic of a TGM325-X Single-Grating Excitation Monochromator




Figure 2-9: DTGM325-X Double-Grating Excitation Monochromator

Swing Mirror Assembly

If three lamps are required on a single-grating monochromator, a swing mirror assembly (as shown in Figure 2-10 below) is attached to entrance port "A", allowing for a third entrance port, entrance port "C".




Figure 2-10: The TGM325-X Swing Mirror Assembly.

2.2.2. Emission Monochromator

The single-grating emission monochromator has one entrance port connected to the sample chamber and two exit ports that can be connected to different detectors.




Figure 2-11: TGM325-M Single-Grating Emission Monochromator

Unlike excitation monochromators, the shutter on the emission monochromator is not computer-controlled. Instead, it is directly interlocked to the sample chamber lid and front hatch. This ensures the detectors are protected against ambient light exposure.




Figure 2-12: Schematic of Single-Grating Emission Monochromator

A computer-controlled swing mirror selects the detector in use. The standard emission monochromator can be fitted with two detectors.




Figure 2-13: DTGM325-M Double Grating Emission Monochromator

A special single-grating emission monochromator is available for imaging applications. This monochromator has exit port B fitted with the standard slit assembly, whereas exit port A has no slit and is reserved for a CCD camera. The swing mirror assembly is in a different position so that the image for port A is outside the monochromator body, at the active area of the CCD camera.

2.2.3. Monochromator Gratings

The monochromator grating determines the spectral coverage and spectral resolution (grating groove density), the efficiency of the optical throughput (grating optimisation wavelength / blaze), stray light and polarisation effects (type of grating – holographic or ruled). The specifications for the standard gratings for excitation and emission monochromator, respectively, are shown in the table below:

Monochromator	Groove Density	Mechanical Wavelength Coverage	Optimisation Wavelength	Dispersion
single-grating excitation monochromator	1800 grooves/mm	200 nm – 900 nm	250 nm (holographic)	1.7 nm/mm
single-grating emission monochromator	1800 grooves/mm	200 nm – 900 nm	500 nm (ruled)	1.7 nm/mm
optional double-grating excitation monochromator	1200 grooves/mm	200 nm – 900 nm	400 nm (ruled)	1.35 nm/mm
optional double-grating emission monochromator	1200 grooves/mm	200 nm – 900 nm	400 nm (ruled)	1.35 nm/mm

Some applications may require additional or alternative gratings. In particular applications in the ultraviolet or in the near infrared spectral range require gratings of different specifications. Please contact Edinburgh Instruments for details on optional gratings.

2.2.4. Slits

Both entrance and exit slits of the two monochromators are continuously adjustable between 0 and 11 mm. With a linear dispersion of 3.5 nm/mm this translates to a spectral bandpass of up to 30 nm:

$$\text{bandpass [nm]} = \text{slit width [mm]} \times \text{dispersion [nm/mm]}$$

For very small slit widths the above formula will be no longer applicable as the monochromator will operate at the resolution limit. At the resolution limit the nominal bandwidth (displayed in the software) may be further reduced. However, this will not further reduce the bandpass, instead it will only reduce the amount light passing through the system.

2.2.5. Shutters

Both monochromators have built-in shutters. The shutters in the excitation monochromator are controlled by the software: The sample is exposed to the excitation light only when measurements are taken or the signal intensity is viewed by the user. At all other times the shutter is closed. This minimises photobleaching and photodegradation of the sample.

The emission shutter is interlocked to the sample chamber lid. The shutter will block when the lid is open or when the sample module is removed from the instrument. This way the detector is protected against potential over-exposure with ambient light.

2.2.6. Filter Wheel

Both excitation and emission monochromators are fitted with filter wheels. When the wheels are in use the software automatically selects the correct filter to remove potential unwanted higher orders. The nominal wavelength that is selected in the software is only strictly correct if the automated filters are in use. If they are not in use higher diffraction orders will pass through the system. These will be detected e.g. at half the wavelength of the excitation, causing scattered or unwanted sample excitation, often leading to misinterpretation of spectra.

The operation of the filter wheels may be disabled from within the spectrometer operating software.

2.3. Sample Chamber

The central component of the FLS1000 Series spectrometer is the sample chamber. This is a large hexagonal box sub-divided into the upper sample chamber, and a lower section containing control electronics.

The sample chamber has two access ports for routine access to the sample and mechanical/optical components inside the compartment: A circular lid that covers a 230 mm diameter round access hole on the top of the chamber and a 260 mm wide front hatch. Both access ports are interlocked, i.e. as soon as one of the two lids is lifted the shutter in the emission monochromator will close to protect the detector(s) against ambient light. Moreover, the lid atop of the sample chamber can be opened and rotated thanks to a hinge at the back of the instrument.

Four slots on top of the sample chamber accept filter holders for excitation beams and emission light. They are not interlocked as the apertures are small and the ambient light entering the sample compartment is minimal.

To minimise stray light, the interior of the sample chamber is coated with a highly absorbing optical grade black paint. Take care not to damage the fragile paint or spill liquids or powders into the sample chamber as this will be difficult to clean off without adversely affecting the paint.




Figure 2-14: Sample Chamber, shown with Optional Polariser and 2nd Attenuator at the Laser Port

The sample chamber has four beam access ports, each 90° apart. This allows for a variety of configuration options such as "T"- geometry (with two emission arms pointing in opposite directions) or "X"-geometry with two excitation arms in addition to the two emission arms. The optical beam height within the sample compartment is 82 mm.

With the sample holder removed, access to the sample compartment is also given from underneath the spectrometer. This enables easy access for the cables or tubes needed to operate the spectrometer with various cryostats.

Attenuator

The standard sample chamber contains an attenuator. This is a circular, continuously variable, reflective neutral density filter.

Reference Detector

The standard sample chamber contains a reference detector behind the attenuator. The reference detector monitors a fraction of the light incident at the sample position. The output of the reference detector is used to correct excitation spectra.

The beam splitter that diverts a fraction of the excitation beam to the reference detector is specially mounted in a way that the angles in respect to the vertical polarisation plane and to the horizontal beam polarisation plane are identical. This eliminates polarisation effects that would otherwise interfere with the correction.

Sample Chamber Options

The sample chamber in Figure 2-14 shows the standard sample chamber, fitted with standard components as well as a popular option of *polarisers* at the excitation and emission beam path and 2nd attenuator at the laser port.

A special version of the sample chamber with a larger beam access port at the excitation side is available. This is needed to mount a vacuum UV refocusing flange for spectroscopy with VUV excitation.

Accessing the Optical Compartments inside the Sample Chamber

Usually there is no need to access the optical components inside the large hexagonal sample chamber. However you may sometimes want to rearrange components for special experiments.

Procedure 1: Accessing the optical compartments inside the sample chamber

- 1-1 Remove the round lid, front hatch and filter holders.
- 1-2 Unscrew the four screws on the top of the large sample chamber cover plate.
- 1-3 Lift the lid carefully rotating it backward. The lid will sit on hinge mechanism.

2.3.1. Standard Lens Optics

The standard sample chamber has lens optics for focussing the excitation light onto the sample and collection of the emission light (see Figure 2-14). The lens assembly has two plano-convex fused silica lenses. In normal operation the lens surface facing the sample is 70 mm away from the centre of the sample position. The demagnification of the lens assemblies is 2.3. A typical slit opening of the excitation monochromator of 0.55 mm (~1 nm bandpass for a standard grating) would produce a 0.24 mm wide beam at the sample position. The typical height of the excitation spot is 3 mm – 4 mm.

The position of the lenses in the sample chamber are adjustable, although for standard applications readjustment should not be necessary.

Before the lens assembly can be moved towards or further away from the sample it is necessary to loosen the grub screw on top of the lens holder. It is then possible to screw the lens assembly into or out of the mount.

2.3.2. Mirror Optics and Other Beam Steering Options

Some applications require special beam focussing and beam steering options.

The table below gives a summary of those applications. Note that the listed applications also require other options and accessories (beyond the mirrors or prisms in the sample chamber), such as lamps, lasers, special detectors.

Application	Suitable Optics	Where?
VUV excitation (<200 nm)	beam re-focussing vacuum flange	excitation
use of external lasers for sample excitation	beam steering prism	excitation
NIR detection (>1700 nm)	mirror optics or lens optics with special lens material	emission
NIR detection (>2500 nm)	mirror optics	emission
opaque samples in front face geometry	mirror optics	excitation and emission

Please contact Edinburgh Instruments for more details about these.

2.4. Sample Holders

2.4.1. Standard Single Cuvette Holder

The standard sample holder is a single cuvette holder as shown below.




Figure 2-15: Single Cuvette Sample Holder with Lens Assembly, Light Cones and Filters

This holder is suitable for liquid samples measured in standard right angle geometry. The sample holder has a socket for quick exchange with other optional sample holders, a thermally isolated base, a sample holder body with feed-through pipes for liquid cooling and a temperature sensor. A spring loaded roller ensures exact positioning of the sample and avoids scratching of the cuvette. The socket of the sample holder is mounted to the base by four M3 screws recessed in slotted holes, providing some freedom for adjustment along the optical axis of excitation. The socket also has slots to hold filters.

The sample holder can be rotated and secured in multiples of 45 degrees.

This sample holder can be removed and replaced at any time without switching off the spectrometer, "hot-swapped". The cable for the temperature sensor (terminated with a D-type connector) can be removed and refitted to the upper of the two D-type sockets next to the emission channel lens assembly.

The sample holder can be completely removed from the instrument by disconnecting the coolant tubes (quick release fittings), the 9-way connector and (if present) the magnetic stirrer. For a temporary repositioning, simply move the sample holder to the right side inside the sample chamber.

2.4.2. Other Sample Holders and Sample Positioning Options

Other sample holders include 3-position sample turrets, front-facing sample holders for liquids, powders and films, holders with magnetic stirrers and a variety of sample cooling and heating options. Additionally, the sample chamber can accept an integrating sphere, adapters for multi-well plate readers, stopped flow assemblies, titrators, computer-controlled sample positioners and optical fibre accessories. For details on these and other options, please contact Edinburgh Instruments.

2.5. Detectors

The FLS1000 usually includes a high-gain photomultiplier (PMT) detector, suitable for photon counting in both steady state and time-resolved applications.

Photomultipliers operating in photon-counting mode provide the highest sensitivity in the spectral range from 200 nm to 1700 nm; for measurements above around 1700 nm analogue detectors are required.

A variety of photo-detectors can be fitted to the FLS1000 Series. When comparing detectors for sensitivity, two factors need to be considered:

- *Detector Responsivity*, usually published by the detector manufacturers as spectral responsivity of quantum efficiency curves
- *Detector Noise*, determined by the background signal (dark count rate or dark current)

For the sensitivity of the whole spectrometer (as opposed to that of the detector alone) it is important to consider other factors such as the efficiency of the monochromators (grating curves), efficiency of the optics, and power and brightness of the light sources in the spectral range of interest.

For details on detectors, alone and in combination with other spectrometer components, please contact Edinburgh Instruments.

2.5.1. Red Sensitive Photomultiplier in Cooled Housing

The standard photomultiplier, PMT-900, is a side window photomultiplier with a spectral coverage from 200 nm to ~870 nm. This detector comes in a cooled housing that provides an operation temperature for the detector of -20 °C to reduce the dark count rate to a minimum.

The cooled housing comprises the optics for refocusing the radiation from the exit of the monochromator to the photoactive area of the detector, the dynode chain socket, thermoelectrical cooling elements and a fan for heat dissipation. A separate power supply, CO1, provides power and controlling circuits for the cooling elements of the detector head.

The detector bias voltage is supplied by the spectrometer controller. This is typically -1200 V.




Figure 2-16: Detector Assembly, comprising Cooled Head to be connected to the CO1 Power Supply



Do not disconnect the CO1 power supply from the detector when is use, as this can permanently damage the thermo-electrical circuit.

Detector Specification

The spectral response of this detector is specified to extend to 870 nm. At that wavelength the quantum efficiency is reduced about 300 times compared to the maximum at 300 nm.

The dark count rate at the operating temperature of -20°C is typically 100 cps.

The detector's minimum response width for time-resolved (TCSPC) measurements is 600 ps. This is the characteristic width of the detector only. Additional contributions for broadening the temporal profile (such as the pulse width of the exciting light source) may also have an effect on the measured pulse width.




Figure 2-17: Quantum Efficiency of Red-Sensitive Photomultiplier

2.5.2. Other Detectors

Other detectors include extended-range red detectors, microchannel plate photomultipliers (MCP-PMTs), thermo-electrically cooled InGaAs detectors, and liquid nitrogen cooled infrared detectors. For details on these, please contact Edinburgh Instruments.

3. Introduction to Fluoracle® Software

The spectrometer's operating software is Fluoracle. Fluoracle controls the hardware of the instrument and provides functions for data acquisition, data analysis and data presentation.

Fluoracle is protected by a USB HASP key (software protection device –SPD) when installed 'off-line' on a different computer for data analysis or presentation.

The following sections give an introduction to Fluoracle. You may also consult the Fluoracle Online Help (Help > Help Topics).

3.1. Starting Fluoracle

To start Fluoracle, double-click the Fluoracle icon on your desktop or select Start > Programs > Edinburgh Instruments > Fluoracle > Fluoracle.

The Fluoracle splash screen is displayed while the software initialises spectrometer components. This auto-calibration takes approximately one minute, depending on how many computer-controlled components are connected to the spectrometer.

Once initialisation is complete, the Fluoracle main window is displayed with an open Signal Rate window (see Section 3.2.5 *Signal Rate Window* on page 32).

3.2. General Features of the Fluoracle Main Window

The Fluoracle Main Window has three areas: a large *Display Area* for measurement results, a *Menu Bar* and a *Tool Bar* on the top of the screen and a *Status Bar* on the bottom.

3.2.1. Menu Bar

The menu bar offers access to all software features of the Fluoracle software. Note that the list of different dialogues available depends on whether or not a data window is displayed and on whether the data are spectral or lifetime data. The list of dialogues accessed through the menu bar also depends on which sample module is in operation. In addition, some features of the instruments (e.g. use of filter turret, display of chromaticity diagrams and quantum yield wizards) can be disabled in the *Options* menu, and consequently the list of dialogues may be reduced.



3.2.2. Tool Bar

The tool bar has a row of icons for commonly-used commands. Some icons will be greyed-out, depending on the currently selected display area window.

- Scan Icons



Opens the Signal Rate, New Spectral Measurement, New Lifetime Measurement, Re-run Measurement, Batch Measurement windows.

- File Icons



Open, save, import, export or print files.

- Scan View Icons



Change the view style for data windows. Depending on the type of a selected scan, switch the plot between 2D, 3D, text, or color map/contour plot. You can also set plot options such as axes fonts and units, screen colours, and labels.

- Measurement Plot Icons



Zoom in to or out of a plot or show the cursor so that you can select and get data on a specific area. Other icons enable a plot grid, view plot peaks, normalise the graphs (in a contour plot) and view or edit plot properties.

- Measurement Container Icons



Join scans into a single measurement container (multiple scans), split multiple scans or extract individual scans from measurement containers.

- Scale Icons



View the Y-axis in logarithmic or linear format.

3.2.3. Scan Display Area

The Scan Display Area is where Fluoracle scans and windows are displayed. You can have multiple scans open at the same time: use the Window menu to manage how they are displayed.

Right-click on a scan to display a pop-up menu. From here you can view the scan's properties, change the view and set multiple scan options.

Note that the scan display area can get cluttered when many scans are open. To obtain an overview select **Tile Horizontally** or **Tile Vertically** from the Windows menu. A large number of open scan displays can be best handled by minimizing displays that are not in focus.

3.2.4. Status Bar

The Status Bar at the bottom displays information about the current scan and spectrometer settings. What is displayed varies depending on whether a scan is in progress and on your particular component configuration.

Both **Tool Bar** and **Status Bar** can be disabled using the View menu options.

3.2.5. Signal Rate Window

The Signal Rate Window (Setup > Signal Rate) is where you setup the wavelength and other parameters before starting a scan. It is displayed automatically whenever you start Fluoracle and is one of the most referred to windows within the software.




Figure 3-1: Signal Rate Window – principal layout

Property	Description
Excitation Wavelength	Set the excitation wavelength and spectral bandwidth.
Emission Wavelength	Set the emission wavelength and spectral bandwidth.
Intensity	Set optical attenuation of the light incident to the sample.
Source Light path	Select the light source to use for the excitation.
Detector Light path	Select the detector to use.
Detector signals	Displays the current signal rates in counts-per-second (cps).

The Signal Rate Window may have additional options depending on your spectrometer's geometry configuration, polarisers presence, number of monochromators and type of sample holder.

3.3. Scan Types and Scan Properties

Measurement data or analysed data are called **Scans** in the Fluoracle software. A variety of different scan plots are available and can be grouped into either spectral scan plots or time-resolved scan plots. Consequently they could be saved into two different file types, annotated with the file extensions .FS for spectral measurements and .FL for lifetime measurements.

This section of the manual briefly describes the scan plot types and outlines the scan properties that are attached to each scan. The scan properties are available via **File > Properties**, or by right mouse click on the active data window and selection of **Properties**. If the active window contains more than one scan a property container is displayed first and a scan must be selected from the list of scans to view the scan properties.

3.3.1. Spectral Data

Excitation Scans

Spectral scans, measured with fixed emission wavelength and scanning excitation wavelength.





Figure 3-2: Excitation Scan, typical Display and Scan Properties

Emission Scans

Spectral scans, measured with fixed excitation wavelength and scanning emission wavelength.





Figure 3-3: Emission Scan, typical Display and Scan Properties

Synchronous Scans

Spectral scans, measured with excitation and emission monochromators scanning simultaneously, using a fixed offset between excitation and emission wavelengths.





Figure 3-4: Synchronous Scan, typical Display and Scan Properties

Emission Correction Scans

Spectral measurement of the emission beam path: This scan contains the spectral information of the emission monochromator and the detector using calibrated light source(s) and is therefore used to spectrally correct uncorrected emission scans. The step size is typically 0.05 nm and the scan is normalised to 1. Note that only Emission Correction Scans can be used for spectral correction of emission scan data. Conversion of normal Emission Scans into Emission Correction Scans is not possible.




Figure 3-5: Example of an Emission Correction Scan, typical Display and Scan Properties

3.3.2. Time-Resolved Data

Time-resolved measurements are data that were acquired by TCSPC, MCS, Oscilloscope or standard single photon counting (kinetic scan). This scan type is sub-divided into decay measurements and measurements of the instrumental response function (IRF). If the measurement is an IRF then the box "Is Instrument Response" must be ticked (refer to the property box, bottom left: *Is Instrument Response*). Windows containing a single decay measurement can be analysed using *Exponential Tail Fit* analysis; windows containing a single decay and an IRF can be analysed using *Exponential Reconvolution Fit* analysis.

Time Scan




Figure 3-6: Time Scan, typical Display and Scan Properties

Fitted Time Scans and Residual Time Scans

The result of a numerical fit is a set of lifetime parameters, as well as two different types of Time Scans: The fitted function and the residual function. These time scans are typically not shown separately (although they can be) but together with the (raw data) IRF and decay.

The properties of the *Fitted Time Scan* reveal the fit parameters. The properties of the *Residual Time Scan* present additional information, such as the Durbin Watson parameter.




Figure 3-7: Fitted Time Scans (red) and Residual Time Scans (green) with the respective detailed Properties.

Autocorrelation Time Scan type is similar to the Time Scan types outlined above and can be shown in the same graphical format with similar scan property boxes.

3.4. Graphical Presentations and Plot Options

Fluoracle provides 2D and 3D graphics, contour plots for data visualisation, as well as a display of the numerical values for viewing and editing purposes.

2D, 3D, and contour graphics have their own colour, font style and line style sets. These can be set up via the *Plot Options*. The three styles are used for the display of data within the Fluoracle and when the graph is copied into the clipboard for further use in different computer applications (**Ctrl + C**). For 2D graphics, a separate set of colours and line styles may be set up for printing purposes (**File > Print**).

3.4.1. 2D Data Presentation

Nine different 2D profiles can be independently set up in respect to their labels on X- and Y- axes, whether the Y-axis should be scaled from Zero or not and whether the Y-axis is displayed in linear or logarithmic scale as default. The use of different 2D profiles eases the handling of different data types. For example, 2D graphics for spectral scans and kinetic data may be set up with linear Y-axis as default, whereas time-resolved data and plots of fitted decay curves are with a logarithmic Y-axis. It should be noted that in order to properly display data normalized to one in logarithmic scale 'scale from baseline' should be unchecked. This option can be found via **right mouse click > Plot options > Fit-Plots**.

The nine different 2D plot defaults can be set up via **Options > Plot defaults > 2D**. The plot of each individual 2D graph may also be changed **right mouse click > Plot options**.

Figure Figure 3-9 and Figure 3-10 below show typical 2D displays for illustration, a "Spectral" 2D plot and a "Fit Plot", respectively.




Figure 3-8: Setup of default 2D Plot Options




Figure 3-9: Example for 2D Display of Spectral data




Figure 3-10: Example for 2D Display of Time-Resolved data (IRF blue, measured decay red)

Task	Operation
Set up 2D defaults	<i>Options > Plot defaults > 2D...</i>
Modify colours and line styles	<i>Plot options > General</i>
Modify axes labels, scaling from Zero, swap between wavelength/time and channels	<i>Plot options > graph type tab</i>
Switch between linear and logarithmic Y-scale	use  and  from the tool bar
Switch cursor on or off	<i>View > Show cursor</i> , or  . When the cursor is in use the X- and Y-values will also be displayed.
Switch grid on or off	<i>View > Show grid</i> , or 
Show maxima, minima, or both	<i>View > Show peaks</i> , or  . Only available for spectral data.
Modify the sensitivity of the peak search	<i>View > Peak settings...</i> Use this dialogue box also to specify, whether you want to display peaks, troughs, or both.
Zoom	use <i>View > Zoom In / Zoom Out</i> or  / 

3.4.2. 3D Data Presentation

Data display in 3D is available if multiple scans or scan maps are contained in the same active window.




Figure 3-11: Example of a 3D display

Task	Operation
Set up 3D defaults	<i>Options > Plot defaults > 3D...</i>
Change colours	<i>Plot options > Appearance</i>
Change the kind of mesh	<i>Plot options > Mesh</i>
Zoom on X, Y, and Z axes	<i>Plot options > Zoom</i>
Rotate 3D display	left mouse click and drag

3.4.3. Contour Plots/Color Maps

The contour plot and/or color map option is available, if the active display contains more than one spectral scan or time-resolved measurement.




Figure 3-12: Colour Map of Spectral data with three Cross-Hairs and Intensity Colour Scale




Figure 3-13: Contour Plot of Spectral data with three Cross-Hairs and Coloured Contour Lines




Figure 3-14: Contour Plot with filled Colour Map and overlaid (white) Contour Lines




Figure 3-15: Colour Map of time-resolved data with one cross-hair and intensity colour scale

Task	Operation
Set up Contour Graph defaults	<i>Options > Plot defaults > Contour...</i>
Change type of contour graph (see above three examples)	<i>Plot options > Colour Map</i> and <i>Countour</i> checkboxes
Change Colour Map scheme	<i>Plot options > Colour Map</i> > toggle between <i>Jet</i> , <i>Hot</i> and <i>Greyscale</i> colour maps
Change colours, add/remove contour lines, make contour lines invisible	<i>Plot options > Contour Colours</i>
Change resolution on wavelength axes (spectral data)	<i>Plot options</i> > use a slider to select between <i>LOW</i> an <i>HIGH</i> resolution
Change resolution on time axis (lifetime data)	<i>Plot options</i> > use a slider to select between <i>LOW</i> an <i>HIGH</i> resolution
Change scale on intensity axis	<i>Plot options > Manual Height Scaler</i>
Zoom	use <i>View > Zoom In / Zoom Out</i> or  / 
Show live cross-hair and axes bars	<i>View > Show cursor</i> , or  . X,Y and Z values are shown when cross hair is live.
Freeze cross-hair and display in bars	hold Ctrl and click left mouse click
Release fixed cross-hair	hold Ctrl, then click left mouse click on cross-hair
Move cross-hair along one axis	hold Shift
Normalise data in display bars	click  in tool bar

3.4.4. View Numerical Data

This facility enables you to view the numerical values of your scan. You can also edit individual data points.

EmMap_mixture_2(Excitation Emission Map)
Experimental Conditions

Scan : Emission Scan : Ex1 290.00nm(ExEmMap)
Wavelength: 310.00 to 574.00 step 2.00 nm

Ex WaveLength: 290.00 nm
Em Slit Size: 2.00 mm
Em Polariser: No polariser
Dwell Time: 0.500 s

Ex Slit Size: 2.00 mm
Ex Polariser: No polariser
Temperature:

Data / 10e5 (Counts)

nm	0	2	4	6	8	10	12	14	16	18
310	1.778	1.694	1.746	2.096	1.968	1.682	2.840	2.903	2.235	2.238
330	2.171	2.120	1.926	1.929	2.023	1.733	1.484	1.399	1.276	1.157
350	1.006	0.867	0.780	0.681	0.583	0.512	0.451	0.393	0.346	0.290
370	0.245	0.217	0.212	0.223	0.204	0.164	0.142	0.114	0.092	0.075
390	0.067	0.063	0.072	0.087	0.098	0.090	0.072	0.056	0.045	0.037
410	0.029	0.025	0.023	0.025	0.031	0.036	0.038	0.039	0.038	0.044
430	0.054	0.071	0.122	0.197	0.227	0.193	0.159	0.150	0.139	0.108
450	0.090	0.076	0.068	0.065	0.066	0.073	0.101	0.121	0.132	0.128
470	0.120	0.112	0.098	0.081	0.069	0.062	0.050	0.041	0.038	0.034
490	0.032	0.033	0.038	0.040	0.040	0.039	0.039	0.035	0.032	0.030
510	0.024	0.024	0.020	0.017	0.016	0.014	0.013	0.011	0.011	0.010
530	0.010	0.011	0.009	0.009	0.010	0.008	0.009	0.006	0.007	0.006

Figure 3-16: Display of numerical values of a scan

3.5. Data Operation with Fluoracle

Fluoracle includes a suite of tools for operating and manipulating your data, accessed from the Data menu.

The Data menu is displayed when the scan display area has at least one scan window open. Options available depend on whether the currently selected scan window contains spectral or lifetime measurements.

Option	Spectral Data	Lifetime Data	Description
Combine	✓	✓	Combines data by adding, subtracting, multiplying, dividing, or (for spectral data) appending.
Scale	✓	✓	Multiplies measurements by a specified factor.
Normalise	✓	✓	Scales measurements to a specified value
Subtract Baseline	✓	✓	Subtracts a stationary background.
Crop Range	✓	✓	Crop measurements down to a specified range.
Smooth	✓		Smooths spectral data using a binomial technique.
Differentiate	✓		Differentiates spectral data using first, second, or higher order derivatives.
Integrate	✓		Integrates spectral data, from the shortest wavelength to the longest.

For more details, see Section 8, Data Operations and Analysis , on page 90.

3.6. Analysis with Fluoracle

Fluoracle includes a suite of wizards for analysing your data. This is accessed from the **Analysis** menu.

The Analysis menu is displayed when the scan display area has a scan window open. Different options will be available depending on whether the open, active scan window contains spectral or lifetime measurements. Options for spectral scans could be Enabled/Disabled in **Options > Analysis** menu. The table below shows only options available for standard system configuration and standard sample holders. Additional options are available for Integrating sphere accessory, Polariser accessories and Transmission detector (See *FLS1000 Series Reference Guide – Integrating Sphere*, *FLS1000 Series Reference Guide – Anisotropy measurements*, *FLS1000 Series Reference Guide – Transmission measurements*).

Option	Spectral Data	Lifetime Data	Description
Chromaticity Plots	✓		Display a chromaticity plot according to either CIE 1931 standard or CIE 1976 standard
Exp. Reconvolution Fit		✓	Reconvolutes a sample decay measurement based on an IRF, specified fitting range, and specified lifetime parameters.
Exp. Tail Fit		✓	Used to analyse scans containing a single decay measurement.
Autocorrelation		✓	Calculate an autocorrelation (residuals) curve.
TRES Data Slicing		✓	Slice TRES data based on a specified start time, stop time and number of slices.

3.7. Software and Hardware Options

The Fluoracle Options menu lets you configure various aspects of the software and spectrometer components.

Option	Description
Correction Files	Set the spectral correction files to be used for correction of spectral scans. See Section.
Windows Options	Set various options for what to do when opening, appending, or closing scan files.
Measurement Options	Set up whether to allow temperature map creation, and the temperature tolerance value for automatic map creation, display reference data, and filter wheel options.
Analysis Options	Enable background subtraction (lifetime data) and whether to enable integrating sphere and absorption calculations for spectral data.
Signal Rate Display Options	Enable display signals of reference detector and transmission detector accessory (if present)
Plot Defaults	Set plot and plot axis options for 2D, 3D, and contour plots.
Sample Holder Options	Select the type of sample holder to use.

Most of the Options dialogues will be explained in more detail in section 4.2.

For more details on the Plot options refer to section 3.4.

4. Starting Up, Setting Measurement and Analysis Preferences and Shutting Down

4.1. Start-up Procedure

Procedure 2: Starting the spectrometer

- 2-1 Switch on the PH1 unit which supplies power to the spectrometer and the data acquisition electronics.
- 2-2 Turn on the CO1 power supply that comes with the cooled detector.
- 2-3 Turn on the Xe2 xenon lamp. If the lamp is in stand-by mode it is sufficient to push the silver button next to the display to ignite the lamp. If the lamp is fully switched off, turn it on using the mains switch on the back and wait until the display shows "Ready to Start".
- 2-4 The microsecond flashlamp, $\mu F2$, has a mains switch on its back panel. Turn it on if you intend to use the lamp. At this point the lamp will not flash. The flashes will be controlled via the FLS1000 software.
- 2-5 Turn a key at the back of any installed EPL/EPLED. When the EPL/EPLED is ready for operation the "ready" light will be constantly on. The laser operation can be started with the "ready" indicator still flashing; however, at this time the laser will not have been fully stabilised.
- 2-6 Switch on the spectrometer computer.
- 2-7 Double click the Fluoracle software. The software will run through an initialisation process checking the various components of the spectrometer.
- 2-8 The microsecond flashlamp's repetition rate is set up within the FLS1000 software via the F2 set up in **Setup > μF lamp setup** or Ctrl+Alt+U.

Note that all settings can be found when the dialogue boxes are opened, such as wavelengths and band-widths, light path selection, formatting and colour selections of 2D, 3D, contour plot displays, etc are all inherited from your last session.

You will find default settings when you operate the software for the first time.

4.2. Setting up Options and Preferences

When working with the spectrometer for the first time, after changing to a different sample holder, or just to optimise certain parameters to your preferences, you may want to adjust some of the options.

4.2.1. Allocate Correction File(s) to the Light Path(s)

True spectra require accurate spectral correction. The files that are used for the correction were measured during manufacture of your spectrometer and were stored in C:\ProgramData\Edinburgh Instruments\Fluorac\Correction. However, correction files can be re-measured, and they can be stored anywhere on your computer. The allocation to the associated light path is made in **Options > Correction Files** (Figure 4-1). To re-allocate the correction files, open up this dialogue box, highlight the light path and select **Modify...**




Figure 4-1: Dialogue Box "Correction Files"

4.2.2. Select Suitable Window Options

A number of Window options is available (Figure 4-2). During your measurement and analysis sessions you will create many data windows. You can reduce the number of windows created by selecting appropriate options in this dialogue box.




Figure 4-2: Dialogue Box "Window Options"

4.2.3. Selecting Measurement Options

The dialogue Measurement Options sets options related to automated temperature map measurements (using a sample temperature controller or a cryostat), enables the automated use of the order selecting filters in both the excitation and the emission path, and enables the display of reference data.




Figure 4-3: Dialogue Box “Measurement Options”

The **Temperature** option in the upper section is used to set a tolerance level for the automated recognition of temperature maps. When measuring temperature maps the sample temperature is recorded at the end of each individual measurement. The recorded value might not be exactly the set value. If one of the recorded temperatures is outside the expected tolerance level the whole data container will not be recognised as a Map. This will have consequences for some data presentations (the Temperature axis label will not be shown) and data analysis options.

The **filter wheel control option** in the middle section of the dialogue enables or disables the automated used of the filter wheels. Note that the correction file needs to be reallocated if the filter wheel is enabled or disabled for the relevant light path. Typically the user will operate the spectrometer with both excitation and emission filters enabled. However, for very special cases, for instance with measurements of high background and very small signal amplitude, the filters can be disabled. However, be aware that the correction files apply only for measurements taken with filters.

The **Enable Anisotropy measurements** option enables anisotropy measurements dialogues. This option is selected by default when computerised polarisers are present and could be selected for anisotropy measurements by using non-standard polarisers such as dichroic sheet polarisers (Refer to *FLS1000 Series Reference Guide – Anisotropy*).

The lower section of the dialogue box contains the tick-box for enabling the **Display Additional windows**. When the option is activated, additional windows will open during spectral scans, displaying the data of the reference detector and raw data of the emission detector.

4.2.4. Selecting signals to be displayed on Signal Rate dialogue

While it is recommended to have signal rates of both detectors, reference and emission, shown on the Signal Rate dialogue, you might disable the signal rate bar of the reference detector and/or transmission detector (if present).




Figure 4-4: Options for Signal Rate display

4.2.5. Enabling Special Analysis Options

Use this dialogue to enable or disable the analysis options that apply to your system. Some analysis features only apply if the correct sample holder hardware is available. For instance, absolute quantum yield calculations will only be needed if an integrating sphere is available (refer to *FLS1000 Series Reference Guide – Integrating Sphere*). The option for polarisation measurements analysis should be selected for analysis anisotropy measurements performed using non-standard polarisers such as dichroic sheet polarisers (Refer to *FLS1000 Series Reference Guide – Anisotropy*).




Figure 4-5: Dialogue Box “Analysis Options”

4.2.6. Setting Options for Data Export

Measurement data are normally stored in a special binary format and have the extensions “.fs” and “.fl” for spectral measurements and lifetime measurements, respectively. Other options include saving numerical results in ASCII codes or CSV files. Plots can also be exported to the Windows Metafile (WMF)

format. You can export data individually via *File->Export* or by pressing the  button on the toolbar. Using the *Options->Export to ASCII* dialog, you can select export options that will allow you to save data simultaneously in different formats when you click the **Save** button.. Here you can specify data delimiters and the code for end of the line, and also indicate whether to include a list of properties (listed before the measurement data) or not.




Figure 4-6: Dialogue Box "Export to ASCII"

4.2.7. Selecting the Sample Holder

This dialogue is important for all computer-controlled sample holders or sample holders with temperature monitor.

You can electrically connect (or disconnect) one of the sample holders that you own, even without leaving the Fluoracle software and without extra software re-configurations. Your sample holder should have been configured by an engineer or at Edinburgh Instruments, therefore it will be offered as an option in the **Sample Holder Configurations** dialogue. Figure 4-7 shows an example of a system that has a standard sample holder (with **Temperature Monitor**), a thermoelectrically cooled single cuvette holder (**Single Cell TEC Sample Holder**) and a sample holder with 3 cuvette positions (**3 Cell Sample Holder**). The selection is set to the sample holder that is connected as shown in Figure 4-7.

For computer-controlled sample holders you can check if the Fluoracle software and the sample holder module are communicating correctly. By pushing the **Apply** button, an initialisation of the module will take place, at the end of which the dialogue box will remain open. Alternatively, by pushing the **OK** button, the dialogue will close while the initialisation will take place.

When the initialisation has been successful some of the displays and functions in the **Signal Rate** screen, in the **Setup** options and in the **Measurement** options will change in accordance with the features that the new sample holder can offer.




Figure 4-7: Dialogue Box "Sample Modules"

4.3. Shutting Down the Spectrometer

The shutting down procedure is simple. Always follow the main rule: Quit software first, then turn off all hardware.

Procedure 3: *Shutting down the spectrometer*

- 3-1 Shut down the software. Remember next time you start your software you will find all settings that you leave the software with.
- 3-2 Turn off the hardware, typically all linked to the power hub PH1. (As the xenon lamp Xe2 is not linked to the PH1 this must be turned off separately – push the silver button near the display of the xenon lamp. This will stop the light emission, but the fans will remain operating for another 60 s to further cool down the bulb and the power electronics. After this time the lamp is in stand-by and can be re-ignited (at the start of your next working session) instantaneously. You can turn the lamp completely off by using the switch near the mains inlet. Do this after the fans have stopped.
- 3-3 You may wish to leave the photomultiplier cooler on to save time at the next start-up.

5. Making Spectral Measurements

FLS1000-s, FLS1000-st, FLS1000-sm, and FLS1000-stm spectrometers can measure steady state photoluminescence spectra in the ultraviolet to the near-infrared spectral range with single photon counting sensitivity. Other instruments of the FLS1000 series can measure spectra too, but overall sensitivity/performance is somewhat reduced due to the lack of the appropriate excitation source.

Typical scans are excitation or emission scans, where the wavelength is the scanning parameter. Other possible scans are scans for producing correction files, maps, and multiple scans. The table below summarises the scan types available for a standard spectrometer. The table also includes spectral anisotropy scans, as polariser upgrades are popular; the acquisition of those spectral anisotropy measurements is included in this paragraph.

Scan Type	Scan Parameter	Second (Map) Parameter
Emission Scan	emission wavelength	excitation wavelength (fixed)
Excitation Scan	excitation wavelength	emission wavelength (fixed)
Synchronous Scan	excitation and emission wavelengths, simultaneously	excitation-emission offset (fixed)
Emission Map	emission wavelength	excitation wavelength, stepwise increased or decreased
Synchronous Map	excitation and emission wavelengths, simultaneously	excitation-emission offset, stepwise increased
Multiple Scans	repeated excitation, emission, or synchronous scans with optional waiting time between repeats	

The recognition of optional sample holders and accessories by the Fluoracle software will automatically increase the type of scans available. For instance, temperature maps will be possible if cryostats or Peltier-cooled sample holders are detected, spatial (sample) maps will be possible with sample X-Y stages or well plate reader attachment, and maps of concentrations will be available with the optional titrator.

The operating instructions in this section assume that the instrument has been started, the spectrometer is in full operating conditions and samples have been prepared. If this is not the case, refer to section 4.

5.1. Preparing for a Measurement – the Signal Rate Dialogue

The following rule applies to all spectral scans:

- 1) Use the **Signal Rate** dialogue box  (Figure 5-1) to select excitation source (Source Light Path) and emission detector (Em Detector Light Path) and to prepare for a measurement. Note, bandwidth and intensity attenuator can only be setup there and not in any other spectral scan dialogue box.
- 2) It is highly recommended to use the Signal Rate dialogue box to set up measurement parameters, i.e. excitation and emission wavelength for which the signal is expected to be at a maximum. During signal rate display and during recording of spectral scans the signal intensity should not exceed 1.5×10^6 cps. Bigger signals will cause saturation effects and non-linearity. Signal exceeding the recommended level will be shown in red (Figure 5-2).



The detector will not break, even with a reading of 5×10^6 cps (which corresponds to a true photon count rate that is much higher), but it is good practice not to exceed 1.5×10^6 cps.

- 3) The photon count rate can be controlled by means of the attenuator (values between 1 and 100 are possible, note that the scale is not linear) and by means of the spectral bandwidth in both the excitation and the emission beam path. In particular for solid or film samples it is important that the impact of the sample positioning and orientation on the signal rate is to be checked at this point. In these cases the sample position and orientation should be adjusted to obtain the highest possible signal rate.
- 4) Use the scan type dialogue box (quick access via ) to set up or modify scan parameters and to start the scan.




Figure 5-1: Signal Rate window for setting up spectral measurements




Figure 5-2: Signal Rate window with emission signal exceeding the recommended level

During spectral acquisitions the status bar on the bottom of the screen informs about the running activities, such as the number of the current scan, the sample temperature, light path, spectral bandwidth, angle of polarisation, filter in use etc.

Active spectral scans can be aborted by either pushing the *Esc* button or by clicking the  tool button on the top right of the Fluoracle screen. Note that the *Esc* button will only be active if the window containing the scan is in focus, not if a different data window has been highlighted.

Details of the parameter settings of a completed scan can be viewed in the *Spectral Scan Properties*. The properties can be accessed via *right mouse click > Properties* if the active window contains only one measurement or via *right mouse click > Properties > Measurement List* if more than one scan is present in the active data window.

Section 5-9 on page 63 provides best-practice advice on making measurements.

5.2. Emission Scans

Scans with a fixed excitation wavelength and variable emission wavelength are called *Emission Scans*. The emission scan range is typically at a longer wavelength than the fixed excitation.

Procedure 4: To set-up and perform an emission scan...

- 4-1 Use the Signal Rate dialogue box to set up source and emission detector light path and other measurement parameters, i.e. excitation and emission bandwidth and level of signal attenuation.
- 4-2 Close the Signal Rate dialogue box and open the dialogue screen Emission Scan Setup.... Upon opening this dialogue box all scan parameters of the last emission scan will be shown as default, apart from those that have been changed in the Signal Rate dialogue since the last scan.




Figure 5-3: Emission Scan Setup dialogue box

- Change the *Excitation Scan Parameters* in the lower section of the dialogue box.
 - The tab *Correction* offers correction options for this spectral scan. Typically all correction options will be ticked, so that the resulting measurement is the true spectrum, undistorted by instrumental effects. It is however strongly recommended that you measure at some point your sample fully uncorrected, so that you can see impact of the correction on the data.
 - The function *Subtract Background* is switched on by default as soon as one of the automated correction methods has been selected. With the Background subtraction option being selected the background signal will be recorded prior to each scan and before each repeated scan will then be subtracted from the raw scan data..
 - The field *Display Additional Windows* are only present if this option has been chosen in *Options > Measurement Options*. If you wish to record reference detector and/or transmission detector and/or raw (uncorrected) data from emission detector tick the corresponding checkbox. The data will be displayed in separate window/s. This might be useful for subsequent correction.
 - The *Excitation* tab confirms the excitation wavelength and spectral bandpass.
 - The *Emission* tab confirms the spectral bandpass in the emission.
- 4-3 All emission scan parameters can be saved by checking the **Apply** button. This will not start the scan, but it saves the parameter settings should the user decide to close and re-open this dialogue box. By checking the Start button the demonstrated scan parameters will be saved, the dialogue box will close and the spectral scan will start.

- 4-4 The Emission Scan graph is displayed with points plotted in real-time. The auto-scaling Y-axis shows the number of photon counts for the total dwell time measured. If more than one repeat have been chosen, then the data of a new repeat scan will be added to the previous scan. The final scan will therefore display the data for the total dwell time. Note that the count rate display in the status bar at the bottom of the screen is different: This shows the (uncorrected) number of counts per second.
- 4-5 Active spectral scans can be aborted by either pushing the *Esc* button or by clicking the  tool button on the top right of the FLUORACLE screen. Note that the *Esc* button will only be active if the window containing the running scan is in focus.
- 4-6 Once the measurement is finished, you can right-click in the scan window to display a pop-up menu to change view settings and to view the scan properties.
- 4-7 Select *File > Save* to save the scan.
- 4-8 Details of the parameter settings of a completed scan can be viewed in the *Spectral Scan Properties*. The properties can be accessed via *right click > Properties* if the active window contains only one measurement or via *right click > Properties > Measurement List* if more than one scan is present in the active data window.

5.3. Excitation Scans

Excitation scans reveal the absorption properties of a sample. They are useful in particular for the investigation of sample mixtures, as by means of the fixed emission wavelength species can be selectively measured.

Scans where the emission wavelength is fixed and the excitation wavelength is scanned over a pre-defined range are known as *Excitation Scans*.

Procedure 5: To set-up and perform an excitation scan...

- 5-1 Use the Signal Rate dialogue box to set up measurement parameters, i.e. excitation and emission wavelength for which the signal is expected to be at a maximum. Ensure that the correct source (typically the steady state lamp Xe2) and the correct detector have been selected. The photon count rate can be controlled by means of the attenuator (values between 1 and 100 are possible, note that the scale is not linear) and by means of the spectral bandwidth in both the excitation and the emission beam path. In particular for solid or film samples it is important that the impact of the sample positioning and orientation on the signal rate is to be checked at this point.
- 5-2 Close the Signal Rate dialogue box and open the dialogue screen Excitation Scan Setup. Upon opening this dialogue box all scan parameters of the last excitation scan will be shown as default, apart from those that have been changed in the Signal Rate dialogue since the last excitation scan.
- Change the *Excitation Scan Parameters* in the lower section of the dialogue box.
 - The tab *Correction* offers correction options for this spectral scan. Typically the excitation correction option will be ticked, so that the resulting measurement is the true spectrum, undistorted by instrumental effects. As for excitation scans the emission wavelength is a fixed value, the *Emission Correction* will only work like a scaling factor for the entire measurement.
 - The function *Subtract Background* is switched on by default as soon as one of the automated correction methods has been selected. With the Background subtraction option being selected the background signal will be recorded prior to each scan and before each repeated scan.
 - The field *Display Additional Windows* are only present if this option has been chosen in *Options > Measurement Options*. If you wish to record reference detector and/or transmission detector and/or raw (uncorrected) data from emission detector tick the corresponding

checkbox. The data will be displayed in separate window(s). This might be useful for subsequent correction.

- The *Excitation* tab confirms the spectral bandpass in the excitation.
- The *Emission* tab confirms the emission wavelength and spectral bandpass.




Figure 5-4: Excitation Scan Setup dialog box

- 5-3 All excitation scan parameters can be saved by checking the **Apply** button. This will not start the scan, but it saves the parameter settings should the user decide to close and re-open this dialog box. By checking the Start button the displayed scan parameters will be saved, the dialog box will close and the spectral scan will start.
- 5-4 The Excitation Scan graph is displayed with points plotted in real-time. The auto-scaling Y-axis shows the number of photon counts for the total dwell time measured. If more than one repeat has been chosen, then the data of a new repeat scan will be added to the previous scan. The final scan will therefore display the data for the total dwell time. Note that the count rate display in the status bar at the bottom of the screen is different: This shows the (uncorrected) number of counts per second.
- 5-5 Once the measurement is finished, you can right-click in the scan window to display a pop-up menu to change view settings and to view the scan properties.
- 5-6 Select **File > Save** to save the scan.
- 5-7 Details of the parameter settings of a completed scan can be viewed in the **Spectral Scan Properties**. The properties can be accessed via **right click > Properties** if the active window contains only one measurement or via **right click > Properties > Measurement List** if more than one scan is present in the active data window.

5.4. Synchronous Scans

With synchronous scans both the excitation and the emission wavelengths are scanned simultaneously, with the offset between excitation and emission being either Zero or with the emission wavelength at a fixed offset above the excitation wavelength.

Synchronous scans are useful for the characterisation of complex fluorophore mixtures. For Zero-Offset synchronous scans, those species that have a non-zero overlap between the excitation and the emission spectrum will produce characteristic spectral bands. The width of the bands will depend on the width of the spectral overlap, while the amplitude of the bands is proportional to the absorbance and the

fluorescence quantum yield of the species contained in the sample. For Non-Zero-Offset synchronous scans, species with shift between absorption and emission that is equal to the offset will produce the biggest signal.

Procedure 6: To set-up and perform a synchronous scan...

- 6-1 Use the Signal Rate dialogue box to set up measurement parameters, i.e. excitation and emission wavelength for which the signal is expected to be at a maximum.
- 6-2 Close the Signal Rate dialogue box and open the dialogue screen **Synchronous Scan Setup**. Upon opening this dialogue box all scan parameters of the last synchronous scan will be shown as default.




Figure 5-5: Synchronous Scan setup dialogue box

- Change the **Synchronous Scan Parameters** in the lower section of the dialogue box.
- Set **Offset** parameter on the **Emission** tab.



It is important to note that synchronous scans with an offset of zero are only meaningful for samples that have negligible scattering properties. If the sample scatters significantly, an offset different from zero should be chosen with the spectral bandwidth of both excitation and emission smaller than the offset.

- As for excitation and emission scans, the tab **Correction** provides options for the spectral correction for this spectral scan. Typically all correction options will be selected to obtain the true spectrum, undistorted by instrumental artefacts. However, it is recommended to study the effect of spectral correction by disabling one or the other correction option at a convenient time.
- The function **Subtract Background** is switched on by default as soon as one of the automated correction methods has been selected. With the Background subtraction option being selected the background signal will be recorded prior to each scan and before each repeated scan. This will then be subtracted from the raw scan data.
- The field **Display Additional Windows** are only present if this option has been chosen in **Options > Measurement Options**. If you wish to record reference detector and/or transmission detector tick the corresponding checkbox. The data will be displayed in separate window/s. This might be useful for subsequent correction.
- The **Excitation** tab confirms the spectral bandpass in the excitation.

- 6-3 All synchronous scan parameters can be saved by checking the **Apply** button. This will not start the scan, but it saves the parameter settings should the user decide to close and re-open this dialogue box. By checking the Start button the demonstrated scan parameters will be saved, the dialogue box will close and the spectral scan will begin.
- 6-4 The Synchronous Scan graph is displayed with points plotted in real-time. The auto-scaling Y-axis shows the number of photon counts for the total dwell time measured. If more than one repeats have been chosen, then the data of a new repeat scan will be added to the previous scan. The final scan will therefore display the data for the total dwell time. Note that the count rate display in the status bar at the bottom of the screen is different: This shows the (uncorrected) number of counts per second.
- Note that for synchronous scans the parameter on the X-axis is the excitation wavelength. The offset of the emission is registered in the scan properties.
- 6-5 Once the measurement is finished, you can right-click in the scan window to display a pop-up menu to change view settings and to view the scan properties.
- 6-6 Select **File > Save** to save the scan.
- 6-7 Details of the parameter settings of a completed scan can be viewed in the **Spectral Scan Properties**. The properties can be accessed via **right click > Properties** if the active window contains only one measurement or via **right click > Properties > Measurement List** if more than one scan is present in the active data window.

5.6. Emission Maps

Emission Maps are a series of emission scans (see Section 0 on page 55) with excitation wavelengths being systematically changed. Emission maps are often displayed by means of a contour plot graphic; this view is often referred to as EEM (Excitation-Emission-Map)

Procedure 7: To set-up and perform an Emission Map measurement...

Have sufficient background knowledge about the sample's emission properties. It is useful to perform one or more emission scans (refer to section 0 on page 55) with the aim to find the set of excitation and emission wavelengths at which the signal is at the maximum.

- 7-1 Open the **Signal Rate** screen; enter the set of excitation and emission wavelength for the expected maximum signal. It might be appropriate to choose identical bandwidths for both excitation and emission.
- 7-2 Close the signal rate screen and open the **Emission Map Setup** dialogue box.
- Use the lower section of the dialogue box to enter the appropriate excitation range and emission range, select the step sizes, the dwell time and the number of scans. Remember that the map is generated by performing emission scans with sequential increase (or decrease) of the excitation wavelength. For equal resolution on the X- and the Y axes of a future contour plot it would be appropriate to select identical wavelength step sizes for both excitation and emission. The time for a measurement can be estimated by

$$t = \frac{\Delta\lambda_{exc}}{S_{exc}} \cdot \frac{\Delta\lambda_{em}}{S_{em}} \cdot d$$

where $\Delta\lambda_{exc}$ – excitation range; $\Delta\lambda_{em}$ – emission range; S_{exc} – excitation step; S_{em} – emission step; d – dwell time




Figure 5-6: Emission Map Setup

- Select the **Correction** tab and choose the method of spectral correction.
- Select the **Excitation** tab and decide whether the excitation should be stepwise increased – unchecked tick box, or stepwise decreased – checked tick box **Reverse Excitation Order** (refer to Figure 5-7). Use the latter option when there is a risk for noticeable photo-degradation of the sample during the long measurement time. Starting at longer wavelengths (lower excitation energies) reduces the impact of photo-degradation on the measurement.




Figure 5-7: Emission Map Setup – Excitation

- Select the **Emission** tab. If none of the boxes in this tab are ticked, the scanning will cover the full range as specified in the lower section of the dialogue box. However, this potentially means that the emission will scan across the excitation. For scattering samples the cross-scanning is not desirable, in this case tick the upper of the two boxes. If you also want to avoid picking up the second order scatter, tick the lower box, too. Once a box is ticked you have the option to decide how many nanometres away from the actual excitation wavelength the emission scan should start.
- As your FLS1000 has been fitted with automatic filter wheels to eliminate second order effects, there should be no need to tick the second of the two boxes in the upper section. However, if you operate the FLS1000 without the filter wheel (**Options > Measurement Options**), the option **Avoid 2nd order excitation** should be selected, and should then run the scans uncorrected as the correction files contain the filter responses.


7-3 Start the scan: **Start**.

Figure 5-8: Emission Map Setup – Emission

7-4 Spectral map measurements can be aborted by either pushing the **Esc** button or by clicking the  tool button on the top right of the FLUORACLE screen. Note that the **Esc** button will only be active if the window containing the scan is in focus, not if a different data window has been highlighted. When aborting a map measurement the last active scan will be lost, all previous scans of the map will remain and can be further processed.

7-5 Select **File > Save** to save the map.

5.7. Synchronous Maps

Synchronous Maps are a series of synchronous scans (see Section 5.4 on page 58) with the offset between excitation and emission being systematically increased.

Procedure 8: To set-up and perform a Synchronous Map measurement...

Firstly, one must gain sufficient background knowledge of the sample's emission properties. It is useful to perform one or more synchronous scans (refer to section 5.4 on page 58) with the aim to find the set of excitation wavelength and emission offset for which the signal is at the maximum.




Figure 5-9: Synchronous Map Setup

- 8-1 Open the **Signal Rate** screen; enter the set of excitation and emission wavelength for the expected maximum signal. Then select the spectral bandwidths and the attenuator setting so that the detector signal is appropriate, i.e. < 1.5 million cps). It might be appropriate to choose identical bandwidths for both excitation and emission.
- 8-2 Close the signal rate screen and open the **Synchronous Map Setup** dialogue box.
- 8-3 Use the lower section of the dialogue box to enter the appropriate excitation range and the range for the emission offset, select the step sizes, the dwell time and the number of scans. Remember that the map is generated by performing synchronous scans with sequential increase of the emission offset. The overall measurement might take a long time, so keep **Dwell Time** and **Number of Scans** to a minimum. The time for a measurement can be estimated by

$$t = \frac{\Delta\lambda_{exc}}{S_{exc}} \cdot \frac{\Delta\lambda_{offset}}{S_{offset}} \cdot d$$

where $\Delta\lambda_{exc}$ – excitation range; $\Delta\lambda_{offset}$ – offset range; S_{exc} – excitation step; S_{offset} – offset step; d – dwell time

- 8-4 Select the **Correction** tab and choose the method of spectral correction.
- 8-5 Start the map measurement: **Start**.

5.8. Multiple Spectral Measurements

Multiple spectral measurements allow for repeated measurements of either emission scans, excitation scans, transmission scans or synchronous scans. These measurements are ideal to follow slow kinetic processes or to study photobleaching and other long term stability effects of your sample.

Procedure 9: To set-up and perform a multiple spectral measurement:

- 9-1 Open the **Signal Rate** screen and enter the set of excitation and emission wavelengths for the expected maximum signal. Then select the spectral bandwidths and the attenuator setting so

that the detector signal is appropriate, i.e. $< 1.5 \times 10^6$ cps). Take into account that over the length of the full measurement the signal may be growing, or that peak excitation and emission wavelengths may be shifting.

- 9-2 Close the signal rate screen and open the **Multiple** scan setup dialogue box.




Figure 5-10: Multiple Scan Setup

- 9-3 Decide whether you want to produce a series of emission scans, excitation scans, transmission scans or synchronous scans. Once selected by the radio button, you will be able to fine-tune the respective scan parameters through the button **Scan Setup...**. These settings are inherited from the previous single scans, but it is advisable to check them again before committing to the longer multiple scan.

A particular point to consider is the duration for a single scan. You may want to change the dwell time and/or the number of repeats to obtain an acceptable (short) measurement time.

- 9-4 Now select the parameters for the full series of your measurements: how often you wish to repeat the measurement (**Measure the sample ...**) and the time you want to pause between individual measurements (**Wait for ...**)
- 9-5 Start the map measurement: **Start**.
- 9-6 Multiple spectral measurements can be aborted by either pushing the **Esc** button or by clicking the  tool button on the top right of the FLUORACLE screen. Note that the **Esc** button will only be active if the window containing the scan is in focus, not if a different data window has been highlighted. When aborting a multiple spectral measurement the last active scan will be lost, all previous scans of the map will remain and can be further processed.
- 9-7 Select **File > Save** to save the measurements.

Tips for Making Good Spectral Measurements

It is always useful to know your sample's absorption properties before attempting spectral emission measurements. If proper spectroscopic data is not available, just looking at the sample can provide valuable information. For example, samples that are transparent will most likely absorb in the UV, and emit in the blue spectral range. Samples that have yellow colour will absorb in the 400 nm region and emit green or orange, samples that are blue will absorb at around 600 nm – 700 nm and have dark red or even infrared emission.

You should also gain experience with sample concentration and alignment issues. Liquid samples are easier to measure than solid or film samples, as typically less scattering is involved and alignment effects are negligible. However, even with liquids concentration effects can seriously affect the quality of spectral scans. For samples with absorption and emission in the visible spectral range it is useful to physically observe the sample when it is excited:

1. Open the sample chamber lid. This will close the detector shutter (for protection) and will allow you to observe the sample.
2. Select the **Signal Rate** box. This will open the excitation shutter so that excitation light is permitted to the sample.
3. Set the excitation wavelength to a value that is consistent with the sample's main absorption.
4. Look at the sample! Observe the excitation, distinguish between excitation and emission (a filter in front of the eye might help), and gain an insight into potential problems with sample alignment and concentration.

The following is a list of potential effects, problems and sources of errors, together with some tips on how to overcome or minimize them.

Tips for all types of spectral scans:

1 Inner Filter Effect

When working with liquid samples, sample concentrations should be chosen that are small enough to avoid the inner filter effect. If the concentration is high, and a wavelength of excitation is selected that is in a range where the sample has high absorbance, only the cuvette surface facing the excitation beam is strongly emitting and no excitation light might reach the cuvette centre. However, the cuvette centre is imaged into the emission monochromator. Therefore not much light might reach the detector.

It is instructive to study this effect by observing the situation using your eyes, for example using a highly concentrated sample of fluorescein with excitation at about 490 nm. The emission can be observed with the sample chamber lid open and the signal rate dialogue box open (so that the excitation shutter is open).

Tips:

- o Reduce the concentration of the fluorophore(s).
- o Choose an excitation wavelength that is in a range where the sample absorbance is reduced. (This can cause other problems, such as excitation of unwanted sample species or sample contaminations.)
- o Use micro-cuvettes or triangular cuvettes.
- o Use a front face sample holder, this allows to detect the emission at the same surface that also faces the excitation beam.

2. Higher Order Effect

A large spike may appear at a wavelength that has two times the value of the excitation wavelength.

This is the second order of the excitation; the problem is a typical phenomenon of a grating monochromator. In this case the problem is generated in the emission monochromator which, when set to a specific wavelength, will not only transmit this selected wavelength, but also a wavelength with half (and a third, and a fourth...) of the numerical value. The problem also appears in the excitation monochromator, but there it is less frequently seen as the excitation wavelength is often shorter than the emission wavelength and often the second order is highly or even completely suppressed.

Tips:

- Use the integrated filter wheels in your spectrometer. The filter operation is selected (*Options > Measurement Options*). Please ensure you have selected the right correction file.
- Use a wavelength cut-off filter (blocking at shorter, transmitting at longer wavelength) in the emission beam path (behind the sample). The filter's cut-off wavelength should be a little shorter than the starting wavelength of the scan.
- If a scan is to be performed over a wide range (say 300 nm – 800 nm) one might even see the second order of the emission (in this case the emission from 300 nm – 400 nm re-appearing between 600 nm – 800 nm). In this situation the best solution is to select the filter wheel option. A more laborious way is to measure the emission spectrum in two steps, with two different cut-off filters, and subsequent joining (*Data > Combine > Append*).

3. Weak Sample Emission

The sample shows very weak emission. What options are there to increase the signal?

Tips:

- The most important task is to check the positioning of the sample in the sample compartment. A visible inspection of the sample excitation is recommended (look at the sample when the **Signal Rate** screen is open). If the excitation is outside the visible spectrum it might be changed temporarily to 500 nm for an assessment of the situation.
- The sample volume of a liquid sample might be too low. Increase the volume, use a micro-cuvette, or raise the cuvette.
- Bulk or film samples might have to be re-positioned.
- Increase the spectral bandpass. If this is only done in the excitation it is still possible to obtain good spectral resolution for the emission scan. In the worst case both spectral bandpasses in excitation and emission have to be increased.

Please note: With very wide spectral bandwidths there is an increased risk to measure stray light that passes through the monochromators. There is also an increased risk to accidentally over-expose the detector when scanning over the excitation or over the second order of the excitation.

4. Detector Saturation

If the detector is exposed to too much signal (either sample emission, or scattered excitation light), the spectrometer will be unable to “count all photons” that impact on the detector. Although this saturation effect is usually caused by more than one problem, the common term is “detector saturation”.

If detector saturation occurs it will not necessarily harm the detector, but the data will be distorted. It therefore is good practice to stay out of the range of detector saturation.

The biggest cause for “detector saturation” is the randomness of the arrival times of photons at the detector and the subsequent randomness of the electrical pulses produced by the detector. A 100 MHz counter would be able to count 100 million photon pulses equi-spaced in time, but with random arrival times the narrowest spaced pulses may be missed due to the detector dead time, and therefore less than 100 million will be recorded.

Detector saturation is easily overlooked when measuring with short dwell times. For example, a spectral scan that was acquired with a dwell time of 0.01 s will be affected by saturation although the biggest signal is only 15,000 counts. Remember the typical saturation threshold is at 1.5 million counts per second!

Tips:

- Use the Signal Rate screen prior to measurements to check for the maximum possible signal
- Start with narrow spectral bandwidths before attempting bigger bandwidths.
- Use the iris attenuator, if appropriate
- Investigate upper signal limits on scans with different dwell times.

5. Repeat Scanning

The dwell time, multiplied by the number of scans, results in the overall integration time per step (when the scan has finished). Is it better to use longer dwell times or a higher number of scans to improve the signal-to-noise ratio?

Tip:

- Use a shorter dwell time, (e.g. 0.1 s or 0.2 s) and a higher number of scans. This has the advantage that the overall spectrum is revealed at a fraction of the overall measurement time. An unwanted spectrum can therefore be aborted much earlier. Also, slow lamp fluctuations or potential detector background drift will be less significant in the measurement result as a result of multi-scan averaging. It is obvious that setting the dwell time to very short values has limitations. This applies in particular to analogue detectors with built-in lock-in and the measurement of samples with very narrow spectral features. If using these the dwell time should not be shorter than 0.3 s.

6. Spectral Correction

Which type of spectral correction is best? Automated or manual? Which type of automated?

There is no general answer to those questions, it often depends how demanding the specific task is. However, in the majority of applications the full automated correction as described in the sections above is appropriate. Nevertheless, it is unavoidable to study the correction and check out the various options.

Tips:

- Do fully uncorrected measurements and use the manual correction (*Data > Correction*) to correct your data.
- Study the impact of the background (high background, small signal levels).
- Data that are measured in spectral regions of reduced detector sensitivity are significantly modified (in a positive and negative sense) by the correction. Spectral correction in this spectral region is particularly demanding.

7. Samples with strong Absorption-Emission overlap

With many samples the emission spectrum overlaps the excitation spectrum to some degree.

Tip:

- In order to measure the full emission spectrum it might be appropriate to excite the sample not at the peak of the absorbance, but some 10 nm to 50 nm before the peak. This increases the range over which the emission can be scanned without distortions by scattered light.

8. Manual measurements of spectral anisotropy

Anisotropy measurements can be performed manually by using non-standard polarisers such as dichroic sheet polarisers and then performing "normal" scans with the polarisers in the correct orientation.

Tips:

- If polarisers are used that are not recognised by the software, the polariser orientation can be edited manually in the properties of the measurement. Access the file properties (right mouse click on the graph) and select *Properties*. On the bottom right side of the scan properties tick the button *Edit Values*, and then change the polarisation angle from none to either 0° or 90° .
- The option for polarisation measurements analysis should be selected (Refer to section 4.2.5 page 50)
- Once the polarisation angle is edited and provided the active window contains at least I_w and the I_{VH} , the anisotropy can be calculated using *Data > Anisotropy*.

6. Making Time-Resolved Measurements

All spectrometers of the FLS1000 series can acquire time-resolved data, with different time resolution. The time resolution depends on whether the spectrometer can perform standard Single Photon Counting (SPC), Multi-Channel Scaling (MCS) or Time-Correlated Single Photon Counting (TCSPC). All three techniques are single photon counting techniques. This guarantees highest detector sensitivity for all time ranges relevant for photoluminescence phenomena. Combined spectrometers (i.e. FLS1000-st, FLS1000-sm, FLS1000-tm, and FLS1000-stm) use more than one of the above techniques and therefore have an increased range of time resolution. You can find out which of the techniques is applicable to your spectrometer by studying section 1.1 on page 4.

Fluoracle automatically selects one of three techniques to acquire time-resolved measurements. The user only has to select the appropriate excitation light path.

Technique	Time Range						
	pico	nano	micro	mili	sec	min	hour
Standard Single Photon Counting (SPC)				✓	✓	✓	✓
Multi-Channel Scaling (MCS)			✓	✓	✓		
Time-Correlated Single Photon Counting (TCSPC)	✓	✓	✓				

FLS1000 spectrometers can also be upgraded to lifetime data acquisition by means of an oscilloscope (for performing time-resolved measurements with analogue detectors in the near infrared spectral region). This fourth data acquisition method is optional; please contact Edinburgh Instruments for details.

Standard Photon Counting for Kinetic Scans (Milliseconds – Hours)

Kinetic Scans refer to measurements with a fixed excitation and a fixed emission wavelength, while the signal, in photon counts per unit time, is recorded over a time period of typically seconds to hours.

Throughout the measurement the sample excitation may be continuously on, for instance in the study of photo-degradation, or the sample excitation may be on for a selected period of time, for instance in the study of after-glow processes or for the study of both rise and decay kinetics.

The “Kinetic Scan” option is available for all type of instruments of the FLS1000 series. The typical excitation source is the Xe2.

MCS (Microseconds – Seconds)

Multi-Channel Scaling (MCS) is a single photon counting technique that is used for recording luminescence decays in the time range between microseconds and seconds. A pulsed light source with low repetition rate is used to excite the sample.

Standard FLS1000-m, FLS1000-sm, FLS1000-tm and FLS1000-stm spectrometers use the microsecond xenon flashlamp (μ F2) as light source for multi-channel scaling experiments. This lamp has a pulse width of $<2 \mu$ s and a pulse repetition rate between 0.1 Hz – 100 Hz.

For a more detailed description of the MCS technique and alternative excitations sources please contact Edinburgh Instruments.

TCSPC (Picoseconds – Microseconds)

Time-Correlated Single Photon Counting (TCSPC) is a single photon counting technique that is used for recording fluorescence data in the time range from picoseconds to microseconds. A high repetition pulsed light source is used to excite the sample and the photons emitted by the sample are processed using the TCC2 card in the computer.

Standard FLS1000 spectrometers with TCSPC capability use the nanosecond flashlamp (nF920) as the default light source for TCSPC acquisitions. This lamp has a pulse width of <2 ns and a pulse repetition rate of typically 40 kHz.

Picosecond pulsed diode lasers and picosecond pulsed light emitting diodes (EPLs and EPLEDs), mode locked lasers and supercontinuum lasers are also increasingly used as light sources for TCSPC measurements with the range of FLS1000 spectrometers.

For a more detailed description of the TCSPC technique and alternative excitation sources refer to the *FLS1000 Series Reference Guide – TCSPC*.

The following measurement types are available:

Measurement Type	Technique	Parameter	Second (Map) Parameter
Kinetic Scan	SPC	time (milliseconds – hours)	
Time-Resolved Measurement	MCS	time (microseconds–seconds)	
Time-Resolved Measurement	TCSPC	time (picoseconds–microseconds)	
Multiple TR Measurements	MCS	time (microseconds – seconds)	
	TCSPC	time (picoseconds – microseconds)	
Excitation TRES	MCS	time (microseconds – seconds)	excitation wavelength
	TCSPC	time (picoseconds – microseconds)	excitation wavelength
Emission TRES	MCS	time (microseconds – seconds)	emission wavelength
	TCSPC	time (picoseconds–microseconds)	emission wavelength
Time-Resolved Anisotropy	TCSPC	time (picoseconds – microseconds)	measurement of vertical and horizontal polarisation plane

The recognition of optional sample holders and accessories by the Fluoracle software will automatically increase the measurement types available. For instance, temperature maps will be possible if cryostats or Peltier-cooled sample holders are detected, spatial (sample) maps will be possible with sample X-Y stages or well plate reader attachment.

The operating instructions in this section assume that the instrument has been started, the spectrometer is in full operating condition and samples have been prepared. If this is not the case, refer to section 4.

The following sections will instruct how to perform the different types of time-resolved measurements listed in the table above, using the Fluoracle spectrometer operating software. The following rule applies to all time-resolved measurements:

Use the **Signal Rate** dialogue box  to prepare for a lifetime measurement, then use the dialogue box for time-resolved measurements (quick access via ) to set up or modify scan parameters and to start the scan.

It is important to routinely use the **Signal Rate** dialogue box, as this contains some items that can only be setup there and not in any of the dialogue boxes for time-resolved measurements. This applies to bandwidth settings, attenuator settings, polariser settings.

It is essential to pay attention to the signal rate levels, time-resolved measurements typically have lower saturation limits to spectral measurements. These saturation limits depend on the rate of the exciting pulsed excitation source, For TCSPC the signal rate screen therefore contains the information about the repetition rate of the source.




Figure 6-1: Signal Rate window for setting up time-resolved measurements

During the time-resolved measurements the status bar on the bottom of the screen informs about the running activities, such as wavelengths and bandwidths settings, the sample temperature, angle of polarisation, etc.

Active lifetime measurements can be aborted by using the **"Stop"** button of the active acquisition dialogue box.

Details of the parameter settings of a completed measurement can be viewed in the **Time Scan Properties**. The properties can be accessed via **right mouse click > Properties** if the active window contains only one measurement or via **right mouse click > Properties > Measurement List** if more than one scan is present in the active data window.

6.1. Kinetic Measurements

With **Kinetic Measurements** the time course of the signal is followed over a period of time, typically minutes, with a time resolution of a fraction of a second to seconds. Excitation and emission wavelengths remain fixed for the time of the measurement.

The light from the excitation source (which is controlled by the excitation shutter) may be permanently on, or off, or may be switched on and off during the kinetic scan.

Procedure 10: To set-up and perform a kinetic measurement...

- 10-1 Use the **Signal Rate** dialogue box to set up measurement parameters, i.e. excitation and emission wavelength and bandwidth. Ensure that the correct excitation light path and the emission light path have been selected. The photon count rate can be controlled by means of the attenuator (values between 1 and 100 are possible, note that the scale is not linear) and by means of the spectral bandwidth in both excitation and emission. For solid or film samples it is important that the impact of the sample positioning and orientation on the signal rate is to be checked at this point.
- 10-2 Close the **Signal Rate** dialogue box and open the dialogue screen **Kinetic Scan Setup...** Set up the (overall) **Scan Time** and the **Time Resolution**. The combination of both will dictate the total number of data points, which is displayed for reference. This number cannot exceed 10,000.
- 10-3 Set up the options for the excitation shutter. If you select the option **Timed**, you can also decide when the shutter will open, and for how long it will remain open.
- 10-4 Kinetic measurements may be corrected for fluctuations of the excitation light. Select the option **Excitation Correction** if you require this correction. The reference and uncorrected emission detector signals may also be recorded and displayed by ticking the corresponding checkboxes in the field **Show Additional Windows**. This additional option is only available, if activated in **Options > Measurements Options**.




Figure 6-2: Kinetic Scan Setup

- 10-5 Start the measurement: **Start**.

An active kinetic measurement can be stopped by either pushing the *Esc* button or by clicking the  tool button on the top right of the Fluoracle screen. Note that the *Esc* button will only be active if the window containing the scan is in focus.

As with all other time-resolved measurements, kinetic scans can be analysed by exponential fitting. This is available via *Data > Exp. Tail Fit*. However, kinetic processes in the time regime of seconds or minutes often follow have more complex models than standard exponentials, the tail fit option might only provide a “crude” time constant for the overall decay (or rise) process.

6.2. Manual Time-Resolved Measurements – MCS

This section describes how to make a time-resolved measurement in multi-channel scaling mode. This mode requires an excitation source with low pulse repetition. The μ F2 is the standard light source for this. Alternative light sources are available; please contact Edinburgh Instruments for details.

6.2.1. μ F Lamp Setup

The μ F2 microsecond flashlamp is computer-controlled in the way that it is triggered by the spectrometer controller. The unit has an integrated mains power supply that needs to be switched on separately (switch located on the rear panel of the μ F2).



The lamp will only provide optical output (flashes with the repetition rate that has been selected in the μ F Setup window) when measurements are in progress, i.e. when the Signal Rate dialogue box is open or while performing a lifetime measurement (or a spectral scan).

Select **Setup > μ F Lamp Setup** to display the setup window for the microsecond flashlamp.




Figure 6-3: μ F2 lamp setup dialogue box

Property	Description
Lamp Frequency	Set the lamp's frequency from 0.1 to 100 Hz.
Detector Gating	Enable gating (where applicable) to set the trigger delay and gate width.

The pulse repetition rate needs to be matched to the samples decay time. For samples with decay times in the millisecond range a repetition rate less than 100 Hz may have to be selected to avoid background build-up. For samples with decay times in the upper millisecond range, or even with seconds, selection of a very low repetition rate (<1 Hz) is required.

Detector gating circuit is an optional accessory. Please contact Edinburgh Instruments for further details.

6.2.2. Measurement of a Phosphorescence Decay

Photoluminescence decays in the time range of microseconds to seconds are called here (for convenience) phosphorescence decays to distinguish between after-glow or persistent luminescence decays lasting as long as several hours. One of the biggest applications of phosphorescence decays of this time range is the time-resolved emission of lanthanides, and this is actually a special type of fluorescence.

Procedure 11: To set-up and perform a time-resolved phosphorescence measurement...

11-1 Check the signal using the **Signal Rate** dialogue box.

- Ensure that the correct excitation path has been selected (μ F2 or alternative low repetition sources). The operation of a low repetition rate light source can usually be heard.
- Select excitation and emission wavelength and their bandwidth.




Figure 6-4: Signal Rate dialogue of Time-Resolved Measurement by MCS

- Check the detector signal (in cps). Ensure that the signal is not saturated. This is not easy and requires experience, as the saturation limit is not only linked to the lamp repetition rate, but also to the phosphorescence decay time of the sample, which you are about to measure. For a completely unknown sample use no more than 2000 cps – 3000 cps as a signal to start with. Change the amount of signal by adjusting the **Intensity** Setting and/or the bandwidth in excitation and emission.



With the microsecond flashlamp in use, photons will only be generated for the duration of the flashes. Depending on the flash rate of the lamp (typically 100 Hz) and on the lifetime of the sample, photons will only arrive at the detector in 'bunches' and this will result in lower saturation limits.

Start with low signals first, before attempting to go for higher rates!

- Close the Signal Rate dialogue

- 11-2 Open the option for manual lifetime measurements ( > **Manual**). Two dialogue boxes will appear on the right side of the Fluoracle main panel: the upper one provides an extract of Signal Rate options and will be required to start and stop measurements and the lower one is specific for the lifetime settings in multi-channel scaling mode (refer to Figure 6-5).




Figure 6-5: Dialogue box for a manual time-resolved measurement by MCS

11-3 In the upper dialogue box you can start (**New**), Stop, and Pause measurements. You also decide whether the new lifetime measurement is to be shown in a previous data window (**Add** – on the very right side), or in a new data window (**NEW** – very top right). For measurements of the phosphorescence decay the box IRF must be left blank. The lower part in the upper dialogue field is tabbed and some parameters inherited from the signal rate screen can be changed, if needed. However, monitoring of the detector signal as a result of those changes is not possible.

11-4 The lower dialogue box is designed for set up the time-related measurement parameters.

- Time Range** and **Channels** affect the full time window and the number of time bins within the full time range, respectively. Changing one of these two parameters will change the time calibration (in time/channel) that is displayed adjacent to these editable boxes.
- The **Lamp Trigger Delay** parameter will shift the position of the rising edge of the measurement along the time axis. An increase of the trigger delay value will cause a shift to the right.
- The **Stop Conditions** can be chosen between stop after a predefined number after a predefined number of **Sweeps** (light source flashes), stop after a predefined measurement **Time** (in seconds) or stop after a predefined number of **Counts**. If the latter stop condition has been selected, the option will be given to select a certain channel for which this stop condition will be valid. If this option is not chosen (box not ticked), the measurement will be stopped when the peak of the overall measurement has reached the stop condition. Note that with this stop condition option, the measurement may take a little longer as the computer will become busy to search through each data sweep to find the maximum peak count. Select one of the other stop conditions to avoid this potential delay.

11-5 Start the measurement by clicking the **NEW** button. The measurement will stop when either the stop condition is reached, or when the **STOP** button is ticked.

A time-resolved phosphorescence decay measurement is typically analysed by fitting to a model function of one or more exponential parameters. This is available via **Analysis > Exp. Tail Fit**.

Time-resolved phosphorescence decay measurements may also be analysed using numerical reconvolution (**Analysis > Exp. Reconvolution Fit**). However, this is only possible, if the active window contains not only the phosphorescence decay, but also the corresponding instrumental response function (IRF).

6.2.3. Measurement of the Instrument Response Function

In a large number of cases, a measurement of the instrumental response function (IRF) will not be required for the analysis of phosphorescence decay measurements. The width of the IRF will be negligible small, compared to the intrinsic lifetime(s) of the decay, and analysis by the Tail fit routine will be sufficient.

However, for the correct determination of short (<50 μ s) lifetimes, reconvolution will be necessary, and thus the IRF needs to be measured. A scatterer (such as Ludox) may be used instead of the sample. This eliminates the risk that a fraction of photoluminescent light "contaminates" the IRF measurement.

Procedure 12: To set-up and perform a measurement of the IRF...

- 12-1 Open up the **Signal Rate** screen; reduce the attenuator (**Intensity**) to the minimum (1), then change the emission wavelength to a value identical to that of the excitation. This way you will be detecting the (Rayleigh) scattered light from the sample. It is important to note that the attenuator is set to a minimum, and only if no light is detected then it should be stepwise increased while observing the detector signal rate. Reducing the attenuator alone might not even be sufficient to attenuate the signal:



For a correct measurement of the IRF the signal should be no more than **1000 cps** (for a lamp operating at 100 Hz).

- 12-2 Close the signal rate screen and open the Manual lifetime measurement screen. It is not recommended to change any of the parameters after step 12-1 above. Changing sample parameters (such as excitation or emission wavelength) would result in a change of the signal rate that has been critically set up, changing time-related parameters must not be done, as the IRF must be measured under identical time-related conditions as the phosphorescence decay.
- 12-3 Before starting the IRF measurement pay attention that the box **IRF** has been ticked. This way the software will recognise the measurement as an IRF measurement (and not a sample decay). This identification is important for the subsequent data analysis. It is also recommended to add the IRF measurement to the data window that was used for the sample decay measurement, as for the analysis the window must contain the decay and the IRF.
- 12-4 If you forget to tick the box **IRF** before you start to acquire data, then you can still identify the IRF after the measurement has been completed. Access the measurement properties (**right mouse click > Properties > (Properties)**) and tick the box **Is Instrument Response** on the bottom left of the property display.
- 12-5 The IRF is typically measured to the same peak height as that of the decay measurement, but this is not absolutely necessary for the subsequent data analysis. Start the measurement (**New**).
- 12-6 On completion of the IRF measurement this can be used in the analysis of the phosphorescence decay. Refer to chapter 9.3.2 on page 104.

6.3. Manual Time-Resolved Measurements – TCSPC

This section describes how to make a time-resolved measurement in Time-Correlated Single Photon Counting mode. This mode requires an excitation source with high pulse repetition rate. This chapter describes procedures of TCSPC measurements using an EPL light source.

6.3.1. EPL Setup

Select **Setup > TCSPC Lasers Setup** to display the setup window for EPL/EPLED light source.




Figure 6-6: EPL/EPLED setup dialog box for TCSPC

Property	Description
Lasers	Indicates the name of the EPL/EPLED light source
Ex. Wavelength (nm)	Indicates the central wavelength of the EPL/EPLED spectrum.

Select the source you need for lifetime measurements. This parameter has no impact on the operation, it is here for records that will appear in the properties of a measurement

Set the pulse repetition rate of your TCSPC laser (manual switch on the laser head) to 1 MHz. The pulse repetition rate might have to be adjusted later. This depends on the decay time of your sample. With 1 MHz pulse repetition rate most samples with lifetimes from picoseconds to about 20 ns can be measured.

6.3.2. Measurement of a Fluorescence Decay

Photoluminescence decays in the time range of picoseconds and nanoseconds are called here (for convenience) fluorescence decays although this might not always be scientifically correct, as heavily quenched phosphorescence decay processes might also fall into this time range.

Procedure 13: To set-up and perform a time-resolved fluorescence measurement...

13-1 Check the signal using the **Signal Rate** dialogue box.

- a. Ensure that the correct source and detector have been selected.

- b. Ensure emission wavelength has been appropriately selected.




Figure 6-7: Signal Rate dialog for TCSPC measurements

- c. Correct operation of the high repetition rate light source can be verified by the frequency that is displayed on the upper of the two signal rate bars. The pulse repetition rate of the light source should be reviewed taking account of the expected lifetime (or estimating the expected lifetime) of the sample. The pulse repetition period of the light source should not be shorter than 10 times the longest lifetime to be measured.
- d. Check the detector signal (in cps). Ensure that the signal is below saturation and pulse pile-up level. With 1 MHz laser pulse repetition rate the saturation limit (caused by the so-called pulse pile up problem typical for TCSPC measurements) is 50,000 cps (=5% of 1 MHz). Signal exceeding this level will be shown in red. Change the amount of signal by adjusting attenuator on the source and/or the bandwidth in excitation and emission.
- e. Close the signal rate dialogue.
- 13-2 Open the option for manual lifetime measurements ( > **Manual**). Two dialogue boxes will appear on the right side of the Fluoracle main panel: the upper one provides an extract of **Signal Rate** options and will be required to start and stop measurements and the lower one is specific for the lifetime settings in TCSPC mode (refer to Figure 6-8).




Figure 6-8: Dialogue box for a Manual Time-Resolved measurement by TCSPC

- a. In the upper dialogue you can start (**New**), **Stop**, and **Pause** measurements. At the top right you can also select **Add** or **New** whether the new lifetime measurement is to be shown in a previous data window or in a new data window, respectively. For measurements of the fluorescence decay the box IRF must be left blank. The lower part in the upper dialogue field is tabbed and some parameters inherited from the signal rate screen can be changed, if needed. Changes will be made after the **Apply** bottom has been ticked.
- b. The lower dialogue box has a display of the two rates: **Start Rate** = pulse repetition rate from the excitation source; **Stop Rate** = signal count rate from the detector. This signal will be shown in red if it exceed threshold of 5% from Start Rate (Figure 6-9: TCC setup: control over 5% rule for TSCPC measurements).




Figure 6-9: TCC setup: control over 5% rule for TSCPC measurements

Use this dialogue also to set up the time-related measurement parameters:

- c. **Time Range** and **Channels** affect the full time window and the number of time bins within the full time range, respectively. Changing one of these two parameters will change the time calibration (in time/channel) that is displayed adjacent to these editable boxes.
- d. The **Options** box is for an advanced setup of selecting Reverse mode of the electronics, setting up delays, etc. Refer to chapter 6.6 on page 84 for details.
- e. The **Stop Conditions** can be chosen between stop after a predefined measurement **Time** (in seconds) or stop after a predefined number of **Counts**. If the latter stop condition has been selected, the option will be given to select a certain channel for which this stop condition will be valid. If this option is not chosen (box not ticked), the measurement will be stopped when the peak of the overall measurement has reached the stop condition. This is probably the most commonly used stop condition in TCSPC measurements, most measurements will be stopped with the value set to 10000 counts.

13-3 Start the measurement by clicking the **New** button. The measurement will stop when either the stop condition is reached, or when the **Stop** button is ticked.

A time-resolved fluorescence decay measurement may be analysed by fitting to a model function of one or more exponential parameters. This is available via **Data > Exp. Tail Fit**.

Note that the tail fitting procedure is only applicable for fluorescence decays that are not affected by the finite width of the instrumental response function. With a width of the IRF of 0.1 ns - 0.6 ns (detector dependant), tail fitting will only provide satisfactory results for lifetimes greater than approximately 5 ns.

Lifetimes shorter than this must be analysed using numerical reconvolution (**Data > Exp. Reconvolution Fit**). However, this is only possible, if the active window contains not only the fluorescence decay, but also the corresponding IRF.

6.3.3. Measurement of the Instrument Response Function

For the analysis of a fluorescence decay by numerical reconvolution an instrument response function needs to be obtained. The IRF contains the information about the time response of the overall optical and electronic system. Only a correct measurement of this information, and the use of this information in the process of analysis, will result in accurate results and lifetimes that can be trusted.

Procedure 14: To set-up and perform a measurement of the IRF...

- 14-1 Open up the **Signal Rate** screen and change the emission wavelength to a value identical to that of the excitation. This way you will be detecting the (Rayleigh) scattered light from the sample. For a measurement of the IRF the rules of saturation apply, too. The level of signal can be controlled by adjusting the signal attenuator in front of the laser head. For an accurate measurement of the IRF a scatterer (such as Ludox) may be used instead of the sample. This way the risk that a fraction of fluorescent light "contaminates" the IRF measurement is eliminated.
- 14-2 Close the signal rate screen and open the **Manual** lifetime measurement screen. It is not recommended to change any of the parameters after step 12-2. The IRF must be measured under identical conditions compared to the previous measurement of the fluorescence decay.
- 14-3 Before starting the IRF measurement pay attention that the box **IRF** has been ticked. This way the software will recognise the measurement as an IRF measurement (and not a sample decay). This identification is important for the subsequent data analysis. It is also recommended to add the IRF measurement to the data window that was used for the sample decay measurement, as for the analysis the window must contain the decay and the IRF.
- 14-4 If you forget to tick the box **IRF** before you start to acquire data, then you can still identify the IRF after the measurement has been completed. Access the measurement properties (**right mouse click > Properties > (Properties)**) and tick the box **Is Instrument Response** on the bottom left of the property display.
- 14-5 The IRF is typically measured to the same peak height as that of the decay measurement. But this is not absolutely necessary for the subsequent data analysis. Start the measurement (**New**).
- 14-6 Once the IRF measurement it can be used in the analysis of the fluorescence decay. Refer to chapter 9.3.2 on page 104.

6.4. Multiple Time-Resolved Measurements

This measurement option provides a tool to automatically acquire a series of time-resolved measurements (TCSPCS or MCS) by repeating the measurements with or without a preset delay (waiting time). All data is stored together in a single container file.

Typically a Multiple Time-Resolved Measurement is made to establish a possible change in the decay process over a longer period of time. It is recommended that single measurements on this or a similar sample have been made before (refer to section 6.2 on page 74 for MCS measurements and section 6.3 for TCSPC measurements). This way you have a background knowledge on the sample specific settings (such as wavelength and bandwidth) and on time range specific settings (such as the required time range and suitable stop conditions).

Procedure 15: To set-up and perform a Multiple Time-Resolved Measurement...

- 15-1 Open up the dialogue box for multiple time-resolved measurements ( > Multiple). For repeated measurements of an IRF tick the box **Instrument Response**. For a repeated measurement of a sample decay you leave this box un-checked. Define the number of measurements you want to do (**Measure this sample ...**). The box **Save data each time** will remain unticked. This way each of the measurements will start from Zero (With the box checked a new measurement cycle would add the new data to the previous result, effectively continuing the old measurement).




Figure 6-10: Setup of Multiple Time-Resolved Measurements

- 15-2 Before starting the measurement, it is important to review sample and time settings before you start. For Multiple Time-Resolved Measurements the sample specific settings (such as excitation and emission wavelengths and bandwidth, iris settings) are inherited from the Signal Rate window. These settings can be reviewed and changed in the **Sample Set Up** box. In contrast, time specific settings (such as time range, number of channels, stop conditions) are inherited from the Manual measurement dialogue box. These settings can be reviewed and changed in the **Time Set Up** box. A typical oversight is the setup of appropriate stop conditions that applies for the individual measurements.

- 15-3 Now start the multiple measurement: **Start**.

A ongoing multiple measurement can be stopped by either pushing the **Esc** button or by clicking the  tool button on the top right of the Fluoracle screen. Note that the **Esc** button will only be active if the window containing the scan is in focus.

6.5. TRES-Map Measurements

Time-Resolved Emission (or Excitation) Map Measurements is a powerful tool in time-resolved fluorescence spectroscopy. TRES-Map measurements consist of a series of time-resolved measurements (TCSPC or MCS) with the emission (or excitation) wavelength systematically stepwise changed. This Map measurement can be converted to time-resolved emission (or excitation) spectra (TRES). These spectra can provide information that goes far beyond than from (conventional) spectra by steady state fluorescence spectroscopy.

Time-resolved emission measurements are also used for more detailed analysis by numerical curve fitting. For example, systematic changes of the lifetime when changing the wavelength may be revealed. Alternatively, the set of lifetimes may remain constant when scanning over the wavelength, but the amplitude of the lifetime components may systematically change instead. For further information please refer to the advanced data analysis FAST provided by Edinburgh Instruments.

6.5.1. Excitation TRES-Map Measurements

TRES-Map measurements with the excitation wavelength systematically changing are only available when a TCSPC source is used that has a broad spectrum and the wavelength for excitation is controlled by the Fluoracle software. The nanosecond flash lamp nF920 is such a source. Refer to *FLS1000 Series Reference Guide – nF920*.

6.5.2. Emission TRES-Map Measurements

Emission TRES-Map measurements are probably more common than Excitation TRES-Map measurements, as they can be converted into time-resolved emission spectra. In difference to Excitation TRES-Maps they also do not require a broad band excitation source.

Procedure 16: To set-up and perform an Emission TRES Measurement:

- 16-1 Prior to setting up a TRES measurement use the facility for Manual time-resolved fluorescence measurements to perform test measurements (refer to and section 6.3 on page 78 for TCSPC measurements). This way you gain knowledge on the sample specific settings (such as wavelength and bandwidth) and on time range specific settings (such as the required time range and suitable stop conditions) and the time it takes for a single curve to be measured.
- 16-2 Open the dialogue box for setting up TRES measurements ( > *TRES* > *Emission*). Select the range of emission you want to cover and choose the step size. Select a step size that is consistent with the spectral bandwidth in the emission channel. Also, be aware of the number of individual lifetime measurements that is to be performed within the selected spectral range. The excitation wavelength is inherited from the Signal Rate box. Leave the tick-box *Instrument Response* unchecked.




Figure 6-11: Emission TRES-Map setup dialogue box

16-3 Do not start the TRES-Map measurement before you have reviewed the sample and time settings. For the TRES-Map measurement the emission bandwidth is inherited from the **Signal Rate** window. This setting can be reviewed and changed in the **Sample Set Up** box. In contrast, time specific settings (such as time range, number of channels, stop conditions) are inherited from the **Manual** measurement dialogue box. These settings can be reviewed and changed in the **Time Set Up** box.

16-4 Now start the TRES measurement: **Start**.

An ongoing TRES-Map measurement can be stopped by either pushing the **Esc** button or by clicking the  tool button on the top right of the Fluoracle screen. Note that the **Esc** button will only be active if the window containing the scan is in focus.

6.6. Tips for Making Good Lifetime Measurements

Good time-resolved measurements require background knowledge of the sample to be measured, i.e. knowledge about the absorption properties and emission properties. It therefore is desirable that you measure excitation and emission scans first, before attempting lifetime measurements. You are also referred to section 5.9 for tips on making good spectral measurements, starting on page 63 of this user guide.

Time-resolved measurements, however, have specific issues, and they are outlined below:

1. Inner Filter Effect

As for spectral measurements, the inner filter effect should be avoided. This effect will reduce the signal that is available for the lifetime measurement.

If inner filter effect is present, and the sample has a strong overlap between absorption and emission, re-absorption may occur. This can artificially increase the value of the measured lifetime (in the nanosecond time scale).

The same tips apply as outlined in section 8-5 for the inner filter effect.

2. Scatter from the Source

Scattered light from the source might pass the emission monochromator in first or second order, and this interferes with the time-resolved measurement.

This type of interference is recognised by a narrow initial spike “on top” of the decay. With measurements in the microsecond and millisecond regime this narrow initial scatter spike might be just discarded in the data analysis. With decay measurements in the nanosecond time scale and (unwanted) IRF superimposes the fluorescence decay. The region of the IRF might either be discarded in the data analysis, or the additional IRF is numerically treated in the reconvolution process (by adding an additional component with a fixed very short lifetime).

Tips:

- Use a wavelength long wave-pass filter (blocking at shorter, transmitting at longer wavelength) in the emission beam path (behind the sample) to block the scattered spike from the source.
- Discard the “contaminated” region during the data analysis process.

3. Weak Sample Emission

The sample shows very weak emission. What options are there to increase the signal?

Tips:

- The utmost important task is to check the positioning of the sample in the sample compartment. A visible inspection of the sample excitation is recommended (look at the sample when the **Signal Rate** screen is open). If the excitation is outside the visible spectrum, or the light is too weak to be seen by the naked eye, you might change the excitation temporarily to 500 nm for an assessment of the situation. You might even change to another light source (such as the Xe2) for the assessment of the alignment situation.
- The sample volume of a liquid sample might be too low. Increase the volume, use a micro-cuvette, or raise the cuvette.
- Bulk or film samples might have to be re-positioned.
- Increase the spectral bandpass.

Please note: With very wide spectral bandwidths there is an increased risk to measure stray light that passes through the monochromators. There is also an increased risk to accidentally over-expose the detector.

It is worth noting that TCSPC measurements (although they might take very long) can still yield good results even if the signal is as low as the background (a few counts per second). The reason for this is that background photons will be distributed over the full second, whereas the TCSPC data will only be acquired when the time window has been started by the pulse from the source.

4. Pulse Pile-Up on Signal (TCSPC)

Pulse Pile-Up is only a problem in TCSPC measurements, where only one detector count can be accepted for a flash cycle of the excitation source. If more than one detector count arrives at the electronics in a measurement cycle, only the first one is “seen” and the others are ignored. This will

skew the decay measurement towards shorter lifetimes. To avoid this effect the stop rate should never be higher than 5 % of the start rate. For very accurate determinations of lifetimes one should not go beyond 2 %, for “quick and rough” acquisitions 10% might be acceptable.

Tips:

- Attenuate the signal to a level that avoids pulse pile-up (below 5 % of the start rate).
- Use a start repetition rate that is sufficiently high

5. Pulse Pile-Up on Background (TCSPC)

Pulse Pile-Up may also occur on background, if the background count rate is very high, for example with modern NIR-PMTs. The consequence of this type of pulse pile-up is a sloping background level in the resulting TCSPC measurement.

Tips:

- Try to use the detector in a mode that gives the lowest detector background count rate.
- Use a start repetition rate that is sufficiently high.

6. Signal Saturation (MCS)

Multi-channel scaling measurements, particularly those short decays and scatter in short time ranges, can be severely affected saturation. The reason for this is that the number of photon pulses that can “fit into each time bin” is limited. For very narrow time bins this number can be as low as 5 per lamp flash and with the lamp flashing at 100 Hz, the limit for a single bin would be 500 per second.

Tip:

- Measure an unknown sample with no more than 2000 cps – 3000 cps. The signal may then be increased if the sample decays in milliseconds. The signal might even have to be decreased if the sample decay time is only a few microseconds.

7. Advanced option: Measurements in Reverse mode

There are two different acquisition modes in TCSPC measurements performed on FLS1000 system: **Forward** and **Reverse**. They could be set up via **Options** button on the **TCC** dialogue.

In **Forward** Mode the pulse rate from the light source is linked to the START input. This rate (with generally equally spaced pulses in time) is substantially higher than the more or less random pulses from the detector linked to the STOP input (Figure 6-12).




Figure 6-12. TCSPC measurements setup dialogues for Forward Mode

At high rates of the light source the forward mode has a clear disadvantage. Because the vast majority of Time-to-Analog Converter cycles will be started by the START pulse, but never stopped by a STOP signal, it needs to be reset at overflow. The electronics are kept busy more than 20 times the amount actually needed. Therefore Forward mode is recommended for relatively low frequencies of the source i.e. lower than or equal to 1 MHz.

In order to decrease the operating charge of the electronics and to utilize the full capability of signal processing the counting system can be operated in the **Reverse** Mode. In this mode, the signal cable carrying the high-count rate from the light source is linked to the STOP Rate and the low rate is linked to START (Figure 6-13). In reverse mode, the time axis of the memory histogram is internally reversed and the decay is plotted in the same way as for forward mode. Reverse mode is recommended when the source operates at repetition rate from 1 MHz to 100 MHz. Remember while working in Reverse mode you must select a corresponding Time Range i.e. with the name "reverse" on it.




Figure 6-13. TCSPC measurements setup dialogues for Reverse Mode

7. Re-Run and Batch Measurements

7.1. Re-run a spectral measurement

An excitation, emission or synchronous scan, as well as an excitation-emission map or a synchronous map can be repeated automatically (i.e. without manual setting up of the required parameters) using the **Re-run** command. The command is available via the Tool bar () or via right-mouse click > **Re-Run Measurements**.

Note that only spectral measurements performed with standard continuous source (Xe lamp) can be Re-run. The command will be inactive in the following cases:

- If spectral measurement has been performed with any other light source than Xe lamp;
- If spectral measurement has been performed using unrecognised excitation/emission lightpath.

7.2. Batch measurements

Combinations of excitation, emission, synchronous scans, excitation-emission or synchronous maps can be run in **Batch Measurements**. This means that several scans can be set for a sample and measured automatically without the presence of the user. The scans can be set to repeat in loops as many times as required, with a fixed preset delay between each scan. The batch of measurements can be saved and loaded for future use.

Procedure 17: To set-up and perform Batch measurements...

17-1 Open the Batch measurements setup window via **File > Batch Measurements...** menu option or by pressing the Batch icon () on the tool bar.

17-2 The upper panel of Batch Measurements setup dialogue (Figure 7-1) shows all the current content of Fluoracle. The lower panel allows organising a batch protocol for multiple sample measurements. To organise a protocol:

- a. Select a spectral measurement. You can either drag-and-drop it into Batch container or use **Add to Batch** button. Note that only spectral measurements performed with the standard continuous source (Xe lamp) could be added into batch protocol. **Add to Batch** button will be inactive (and drag-and-drop will be not possible) in the following cases:
 - If spectral measurement has been performed with any other than Xe lamp light source;
 - If spectral measurement has been performed using unrecognised excitation/emission lightpath.
- b. Using **Add > Delay** if you wish to set a delay in seconds between consecutive measurements.
- c. To repeat a series of measurements a desired number of times add **Loop (Add > Loop)** before the first one in the series and **End Loop (Add > End Loop)** after the last one.
- d. Change the content and order of measurements within a Batch protocol by using the **Delete, Move Up** and **Move down** buttons.

- e. You can **Save** a Batch protocol into a file of *.fb type. Later you can **Load** a previously saved Batch protocol.




Figure 7-1. Batch Measurements setup dialogue

- 17-3 Press Run to start Batch measurements. The Batch Progress window (Figure 7-2) will appear and the currently performed stage of measurements will be indicated by green triangle. Every measurement will appear in a separate plot window.




Figure 7-2: Batch Measurement progress window.

8. Data Operations and Analysis on Spectral Data

A variety of the data operation tools and data analysis options is available to spectral data. Some data manipulation options are common to both spectral and lifetime data, others are specific. Therefore, the list of options that is available by clicking the pull down menu **Data** will vary depending on the data file type.

A summary of the data operation tools and data analysis options for spectral and time-resolved data is given in section 3.5. "Data Operation with Fluoracle" on page 44 of this User Guide.

This section describes in detail how to use the standard data operation tools for spectral data. A range of non-standard data operation tools is available. They may be activated via **Options > Analysis Options > Spectral Data**. These advanced options require data that can only be measured with special hardware (e.g. integrating sphere, absorption accessory).

8.1. Add, Subtract, Multiply, Divide

With these data options a single spectral scan (of a set of spectral scans) can be added to (subtracted from, multiplied or divided) to another spectral scan.

For example, to add two different curves together, open the dialogue box **Data > Combine > Add...** (refer to Figure 8-1). The list that is provided within this dialogue box contains all data windows that are open or minimised within the Fluoracle program. The sequence is in order of time accessed. In the example of Figure 8-1 the spectral scan **Example Data 2** will be added to the scan **Example Data 1**. Instead of **Example Data 2** another measurement may be picked that may also be in the data container **Example Data Set**. Double click on the container to open it. In order to return to the original list of files and containers, it is necessary click on the  icon that will have opened up on the top right of the dialogue box.




Figure 8-1: Add spectral scans

The spectral scans to be added must have the same step size. It is not necessary that they have the same spectral coverage. If the spectral ranges do not fully overlap the resulting file will extend from the minimum of the two files to be added to the maximum of those files. The data that are not available are automatically treated as Zero.

The resulting curve(s) can be added to the original window, or they can be shown in a new data window. This is selected by the corresponding buttons on the lower part of the dialogue box.

By clicking the **Apply** button the numerical operation will be processed and the new data will be displayed. The dialogue box will remain open for further data adding. Selecting **OK**, in contrast, will perform as with **Apply**, but will automatically close the dialogue box.

8.2. Append

The **Append** operation is only available to spectral data and correction scans. It has been designed for the purpose of joining spectral scans of different spectral ranges into one. The **Append** function will only work if the two different spectral scans to be joined have at least one overlapping wavelength point. The spectral scan that is to be appended on must have a spectral coverage that is at shorter wavelengths than the one that will be appended with.

In contrast to the other **Combine** operations described in sections 8.1, only a single spectral scan can be appended on. If the active data window contains more than one scan the **Append** option will not be available.

Layout and functionality of the **Append** dialogue box are similar to the those described in section 8.1. The difference is additional information and select options that are given beneath the list of data.




Figure 8-2: Append dialogue box

If the selected file fulfils the conditions that are required to append it to the original file the overlap range will be displayed and the option will be given to select a wavelength within the overlapping range (centre of the overlap as default). The selected wavelength in the **Append at** box has no function if the **No matching** option has been selected. However, if either **Scale** or **Offset** have been selected, the data that will be appended to the original curve will be scaled or offset so that they match the original curve at that selected wavelength. The resulting spectrum will have the data of the original curve up to the wavelength selected in **Append at** and will have the scaled or offset data of the second curve from there onwards. An averaging in the overlapping range is not provided.

8.3. Scale, Normalise, Subtract Baseline, Crop Range

Scale

The **Scale** operation multiplies a single measurement or a set of measurements contained in a single data window by a factor.




Figure 8-3: Dialogue box for scaling data

By selecting **Scale** from the **Data** menu the dialogue box as shown in will open up. The scaling factor may be edited and the option for the display of the resulting curve may be selected.

Note that a negative scaling factor is possible; this can be used to invert the data. Note that a scaling factor of Zero is not possible.

This scaling will apply to all curves that are contained in the data window. **Apply** will process and display the data, but will keep the dialogue box open. **OK** will process the data and will automatically shut down the dialogue box.

Normalise

The **Normalise** operation scales active data either to the same peak value (**Normalise > To Peak Value...**, Figure 8-4) or to the same intensity value at a specified wavelength (**Normalise > At Wavelength...**, Figure 8-5). The target intensity value is given in the "to" edit box. This value can be edited manually. Alternatively, the curve is selected from the list of files below. As soon as one measurement is highlighted, the maximum of this measurement (or intensity at a selected wavelength) will be found and will be displayed in the "to" box. If a measurement container is highlighted, the maximum of all curves in that container will be displayed.




Figure 8-4: Normalise Data To Peak Value dialogue box




Figure 8-5: Normalise Data At Wavelength dialog box

Normalisation to negative numbers will invert the data and all data will have a common minimum value. Normalisation to Zero will set all curves to Zero.

Note that when original data are displayed in the logarithmic scale, the resulting data may not be suitable for demonstration in this way, e.g. when curves are normalised to ≤ 1 . The user may have to change the graph to linear scale. In some cases even that may not be sufficient to view the data and the *Scale from Baseline* box might have to be un-checked (*right click on the graph / Plot Options / Spectral / ...*)

Subtract Baseline

This option is for subtracting a stationary background, from a single spectral scan, or from a set of spectral scans.




Figure 8-6: Subtract baseline

The number that is displayed in the editable field will be subtracted from the data of the active data window. Upon opening the dialog box the editable field is normally blank. An exception is a raw single scan that was measured without automatic background subtraction. As the background had been recorded, it is now displayed as a default value for manual background subtraction.

A negative number can also be used in the Baseline value to introduce a higher background in the data.

Crop Range

Spectral scans may be cropped, i.e. the spectral range may be cut back. This function is useful for removing unwanted data from the left or the right side of the graph, in particular if the **Append** function (see section 8.2) is to be used afterwards.

To crop the data range of a single curve or a number of curves in the active data window, open the dialogue box and edit the default values. The default values are the far left and the far right extremes of the current data range.




Figure 8-7: Crop Range dialogue box

8.4. Smooth

Spectral data can be smoothed. The binomial smoothing technique is applied, with the option to include from 1 up to 4 neighbouring data points on either side of each point to be smoothed, corresponding from 3 up to 9 **Number of Points**, respectively.




Figure 8-8: Smooth dialogue box

Note: excessive smoothing may alter real spectral features. Smoothing is best performed on oversampled data, i.e. curves with high number of data points for the spectral features involved.

8.5. Differentiate and Integrate

Differentiate

Use this data operation to differentiate the spectral scans. Derivatives of the first order, the second order, even up to the 4th order can be performed. For the reduction of noise, smoothing can be applied prior to differentiation.




Figure 8-9: Differentiate dialogue box

Integrate

The **Integrate** operation is only available for spectral data. The integration of the whole spectral scan will be performed, starting at the short wavelength side and ending at the data point of the longest wavelength. Thus, the longest wavelength data point of the resulting integrated curve represents the integral number of the whole original curve.

If an integral number is required of a range that starts at longer wavelengths than the first data point of the original scan, then you must either crop the original scan on the short wavelength side before using the **Integrate** operation, or you subtracts the Y-value taken at the long wavelength side of the integrated curve minus the Y-value taken at the short wavelength side.

Note that the **Integrate** and **Differentiate** functions are not fully reversible. Subsequent application of integration and differentiation to a spectral scan will return the original data apart from a scaling factor and a possible Y-axis offset.




Figure 8-10: Integrate dialogue box

8.6. Spectral Correction

This dialogue enables you to perform spectral correction post acquisition, using stored correction files specific to your spectrometer.

Correction files are only available for correction, if they have been correctly assigned for. This assignment is made in *Options > Correction Files....* Maintaining the correct correction file assignment is critical for producing true spectra.




Figure 8-11: Assignment of Correction Files

When the active window contains one or more spectral scans, the option for manual post acquisition correction will be available (*Analysis > Correction...*). The dialogue box shown in Figure 8-12 will open.




Figure 8-12: Dialogue box for Manual Emission Correction

The left panel of the dialogue will show a list of all scans of the active data window. You can select one scan or a group of scans of the same type to correct. To **Select All** tick the checkbox below the list of scans. You should select a corresponding lightpath from the table *Detector Lightpath&Correction files*. Now you can press Ok to perform the correction of your raw data.

8.7. Chromaticity Plots

The chromaticity of a colour is the hue (the property that distinguishes red from purple) and saturation or purity (the property that distinguishes red from pink). It is used to quantify the perception of colour.

The Fluoracle software can produce chromaticity plots according to CIE 1931 and CIE 1976 according to the standards of CIE: Commission Internationale de l'Eclairage – International Commission on Illumination.

The human eye has photoreceptors called cone cells for medium to high brightness colour vision. The Cone cells have sensitivity peaks at blue (420 nm – 440 nm), green (530 nm – 540 nm) and red (560 nm – 580 nm). All perceived colours can be translated into differing amounts of three primary colours corresponding to the cone sensitivities. The different levels of stimulus to the three types of cone cell, or the amount of each primary colour, give three parameters which can describe any perceived colour, these are the tristimulus values, X, Y and Z.

The formulas for the calculation of the tristimulus functions from the measured spectrum and the three standard observer functions are given in the appendix *Formulas and Definitions*, equations 1-2. For the calculation of the CIE 1931 coordinates, x and y, refer to equations 3 and 4. For the calculation of the CIE 1976 coordinates, u' and v' refer to equations 5 and 6.

Procedure 18: To generate a chromaticity plot . . .

- 18-1 Enable chromaticity calculations in analysis options. *Options > Analysis Options > Enable CIE Chromaticity Calculation* (Figure 7-12).




Figure 8-13: Analysis Options with option to enable CIE Chromaticity Calculation

- 18-2 Open the emission scans on which to perform the analysis. Perform the calculations in the analysis menu; *Analysis > Chromaticity*. A warning message will be displayed if the scan(s) are not emission scans.
- 18-3 Scans can be measured with any step size. The colour-matching functions are used to interpolate any further data required.

- 18-4 Using the two tabs on the top of the Chromaticity window, choose between the CIE 1931 and CIE 1976 standard.
- 18-5 You can set reference points. Use the *Delete*, *Edit*, *Add* buttons on the bottom of the screen for this purpose.





Figure 8-14: Chromaticity Plots displayed in CIE 1931 (top) and CIE 1976 (bottom)

9. Data Operations and Analysis on Time-Resolved Data

A variety of the data operations and analysis tools is available to time-resolved data. Some data tools are common to both spectral and lifetime data, others are specific. Therefore, the list of options that is available by clicking the pull down menu **Data** will vary depending on the data type and (for windows containing a set of data) on the properties of the data set.

A summary of the data operation tools and data analysis options for spectral and time-resolved data is given in section 3.5. "Data Operation with Fluoracle" on page 44 of this User Guide.

This section describes in detail how to use the operations and analysis tools that are provided for time-resolved data.

Some analysis routines are numerically complex. For more details on the lifetime data analysis used in the Fluoracle please contact Edinburgh Instruments.

Edinburgh Instruments also offers a software package for advanced lifetime data analysis: *FAST*. This software is sold separately.

9.1. Add, Subtract, Multiply, Divide

These operations are similar to those described for spectral data (refer to section 8.1 on page 90).

However, pay attention that numerical operations like these may have an impact on the noise characteristics of the time-resolved data. For instance, adding one time-resolved measurement to another independently measured time course will maintain the Poissonian noise distribution, whereas adding an identical measurement to the original data will corrupt the Poissonian noise. A correct Poissonian noise distribution is important for fitting time-resolved data.

An **Append** function is not available for time-resolved measurements.

9.2. Scale, Normalise, Subtract Baseline, Crop Range

These operations are similar to those described for spectral data (refer to section 8.3 on page 92).

When using these functions pay attention that this will adversely affect the noise distribution of the time-resolved data. For instance, a non-manipulated time course might be fitted with an excellent fit and a χ^2 value of 1.00. The same measurement, multiplied by a factor 10, can only ever have a minimum χ^2 of 3.16 ($=\sqrt{10}$)!

Another difference to spectral data is the graphical demonstration: Time-resolved measurements are often demonstrated in a semi-logarithmic plot. This plot might not be any more appropriate after the data operation has been applied – the resulting data window might appear to be empty. Change over to linear scale and/or **Scale from Zero** (via **Plot Options**) to retrieve the measurement.

9.3. Analysis of Luminescence Kinetics

Fluoracle software enables to fit raw luminescence kinetics to an exponential components model having up to 4 exponential terms:

$$R(t) = \sum_{i=1}^4 B \exp(-t/\tau_i) \tag{7}$$

The fit will produce lifetimes τ_i , pre-exponential factors B_i , a numerical value for a background level A and . the intensity contribution of each of the exponential components to the overall fluorescence intensity '*Rel,%*'

$$\varphi_i = \frac{B_i \tau_i}{\sum B_i \tau_i} \cdot 100\% \tag{8}$$

The table of results will contain also two additional lifetime values, lifetime averaged by amplitude and lifetime averaged by intensity, that are often used to characterise complex decays. The amplitude averaged lifetime $\langle \tau \rangle_{amp}$ is the lifetime a fluorophore would have if it had the same steady-state fluorescence as the fluorophore with several lifetimes. This value is calculated by:

$$\langle \tau \rangle_{amp} = \frac{\sum B_i \tau_i}{\sum B_i} \tag{9}$$

The intensity average lifetime $\langle \tau \rangle_{int}$ is the average amount of time a fluorophore spends in the excited state. This value is calculated by:

$$\langle \tau \rangle_{int} = \frac{\sum B_i \tau_i^2}{\sum B_i \tau_i} \tag{10}$$

To analyse sample decay measurements that are not (or not significantly) affected by the instrumental response function select an option '*Analysis-> Exp Tail Fit...*'

Alternatively, use the option for reconvolution fit '*Analysis-> Exp Reconvolution Fit...*'

To learn more about the numerical fitting algorithm and about the meaning of the fit parameters, please contact Edinburgh Instruments.

9.3.1. Tail Fit Analysis

Tail fit analysis will only produce meaningful results, if the sample is allowed to decay without further "disturbance" by excitation.

Figure 9-1 below shows – on the example of two different measurements "A" and "B" – the selection of a suitable tail fit region (region with white background). The tail fit region excludes the time period for which the IRF is still present. For obtaining proper fit results, the far end of the tail should be excluded, it only contains detector background and noise and no valuable information on the decay kinetics




Figure 9-1: Tail fit regions, two different examples

Procedure 19: To perform a tail fit...

- 19-1 Ensure the active window contains only one sample decay measurement, or a sample decay measurement with additional IRF.
- 19-2 Inspect your sample decay measurement and decide whether a good model to fit to would be a single exponential decay, or a function that contains more than a single exponential term. Decay processes that are likely to be of single exponential nature are recognised by a sloping straight line in a semi-logarithmic plot (see tail fit region of decay "A" in Figure 9-1). Decays that deviate from a straight line (see tail fit region of decay "B" Figure 9-1) will not be satisfactory fitted with a single exponential term.
- 19-3 Use the **Zoom** function to select the tail fit region.
- 19-4 Select **Data > Exp. Tail Fit**. This will open up a the dialogue box as shown in Figure 9-2. The fitting range shown on top of the dialogue box is inherited from the zoom in step 3 above, provided the zoomed region starts at, or beyond, the peak of the decay. If the zoomed region starts before the peak of the curve, then the **Fitting Range From** value will automatically be set to the peak of the curve. Remember, the peak is not necessarily the most appropriate start of the fitting range. The fit range for a tail fit should start when the IRF is completely finished.
- 19-5 Enter your estimate for the first lifetime, τ_1 . This value will serve as a starting value in the following fitting process. A sufficiently good estimate for τ_1 can be made by looking at the original data: Remember that for a single exponential decay the lifetime τ is the time it takes to decay from 100% to 37% ($0.367=1/e$). Even if your judgment of the decay was that the decay is more complex than single exponential (see step 17-2 above) it is probably a good idea to start with a single exponential fit. The fit result of the single exponential trial will provide you with information about the estimate of a potential second lifetime component.




Figure 9-2: Tail Fit dialogue box

19-6 You do not need to enter an estimation for the background value **A**. If this box is left blank the fit will be started with an initial background of Zero. The background **A** is a numerical value that accounts for a constant detector background or for a very long underlying photoluminescence signal that has built-up as a consequence of the repetition rate of the light source.

19-7 Click **Apply** to start the fit.

19-8 The fit result will be presented almost instantly. The fit result will comprise:

- ✓ an updated fit-plot which contains the raw data, the fitted function (file suffix **F1**), and the residual function (file suffix **F1R**). The plot also contains a display of the fitted lifetime τ_1 and of the fit quality parameter χ^2 .
- ✓ an updated fit dialogue box, containing the calculated lifetime, τ_1 , the calculated pre-exponential factor, B_1 , the asymptotic standard errors for these parameters, the calculated background, **A**, and the goodness of fit (or chi-squared), χ^2 .

19-9 You now need to make an evaluation of the quality of the fitted result, based on two parameters;

- ✓ the value of the χ^2 : For a good fit this should be close to Unity (Refer to chapter 9.4.1 on page 108 for details).
- ✓ the residuals (shown in the fit plot): For a good fit they should be randomly distributed around zero.

Figure 9-3 shows as satisfactory single exponential fit result of decay "A" of Figure 9-1. In contrast, Figure 9-4 (showing a tail fit result of decay "B" of Figure 9-1) is an example for an unsatisfactory fit result.

19-10 Repeat the fitting process, if the fit result was not satisfactory. Add a second exponential lifetime, τ_2 ; this will test for a double exponential model. Repeat with three and four exponential terms, if necessary.

Fit results of models with 2, 3 and 4 exponential components will show a Rel% figure that is different from 100%, provided none of the pre-exponential factors is negative. The Rel% figure states the intensity contribution of each of the exponential components to the overall fluorescence intensity.

19-11 When you are happy with the fit you can print the results directly from within the dialogue box (**Print**), or you can send the data to the clipboard, either as text or as image (**Copy as Text**, **Copy**

as Image), for use in other software applications. To copy the fit plot, close the tail fit dialogue box, then press Ctrl+C on the graph.

Fit results (file suffix *F1*) and residuals (file suffix *F1R*) may be reviewed any time later by looking up the properties of the files with suffices *F1* and *F1R*, respectively. The properties of the residuals file also contain extended fit parameters, such as the Durbin-Watson parameter (Refer to chapter 9.4.4 page 109 for details).




Figure 9-3: satisfactory Tail Fit




Figure 9-4: unsatisfactory Tail Fit

9.3.2. Reconvolution Fit Analysis

In contrast to the Tail Fit Analysis described in the previous section, where only the decay was analysed, Reconvolution Fit Analysis extract the lifetime parameters from the whole time-resolved sample measurement, including the rising edge of the sample response. This way shorter lifetimes, that may even be shorter than the pulse width of excitation, may be reliably extracted.

In order to do numerical reconvolution, not only a measurement of the sample decay is needed, but also an accurate measurement of the IRF.

Figure 9-5 below shows – on the example of two different measurements “A” and “B” – the selection of a suitable fit region suitable for Reconvolution Analysis (region with white background). Note that – in contrast to Tail Fit Analysis – the region of the IRF is included. This way all information of the sample decay can be processed during the analysis. This is clearly beneficial for the analysis of a decay as shown in the example “B”, but even can reveal a short lifetime component of the measurement example “A”.




Figure 9-5: Fit regions for Reconvolution Fit Analysis

Procedure 20: To perform a reconvolution fit...

- 20-1 Ensure the active window contains one sample decay measurement and the corresponding IRF.
- 20-2 Inspect your sample decay measurement and decide whether a good model to fit to would be a single exponential decay, or a function that contains more than a single exponential term. Decay processes that are likely to be of single exponential nature are recognised by a sloping straight line in a semi-logarithmic plot (see decay “A” in Figure 9-1). Decays that deviate from a straight line (see decay “B” Figure 9-1) will not be satisfactorily fitted with a single exponential term.
- 20-3 Use the **Zoom** function to select the tail fit region. This should include the region of the IRF, otherwise refer to the previous section for Tail Fit Analysis.
- 20-4 Select **Data > Exp. Reconvolution Fit...** This will open up a the dialogue box as shown in Figure 9-6. The fitting range shown on top of the dialogue box is inherited from the zoom in step 18-3 above.
- 20-5 Enter your estimate for the first lifetime, τ_1 . This value will serve as a starting value in the following fitting process. A sufficiently good estimate for the start value (for a single exponential fit) is 1/10 of the full scale of the fitting range. You will get experience in finding a suitable start value for τ_1 . It is also good to know that any lifetime estimate that is between 1/10 to 10 times of the fit results will be sufficiently good for a successful single exponential fit. For fits with 2, 3 and 4 exponential terms, however, the estimates have to become more and more accurate.




Figure 9-6: Reconvolution Fit dialog box

20-6 You do not need to enter an estimate for the background value **A**. If this box is left blank the fit will be started with an initial background of Zero. The background **A** is a numerical value that accounts for a constant detector background or for a very long underlying photoluminescence signal that has built-up as a consequence of the repetition rate of the light source.

20-7 Click **Apply** to start the fit.

20-8 The fit result will be presented almost instantly. The fit result will comprise:

- ✓ an updated fit-plot which contains the raw data, the fitted function (file suffix **F1**), and the residual function (file suffix **F1R**). The plot also contains a display of the fitted lifetime τ_1 and of the fit quality parameter χ^2 .
- ✓ an updated fit dialogue box, containing the calculated lifetime, τ_1 , the calculated pre-exponential factor, B_1 , the asymptotic standard errors for these two parameters, the calculated background, **A**, shift of the decay relative to IRF, δ , and the goodness of fit (or chi-squared), χ^2 .

20-9 You now need to make an evaluation of the quality of the fitted result, based on two parameters;

- ✓ the value of the χ^2 : For a good fit this should be close to Unity. (Refer to chapter 9.4.1 on page 108 for details).
- ✓ the residuals (shown in the fit plot): For a good fit they should be randomly distributed around Zero.

Figure 9-7 shows unsatisfactory single exponential fit results for decay "A" and "B", respectively. It is clear from the look of the residuals that a fit with more than a single exponential model function would give better results. Figure 9-8 shows a satisfactory fit result.

20-10 If you are not satisfied with a single exponential fit, repeat the fitting process with an added second, third, or even fourth exponential lifetime.

Figure 9-9 and Figure 9-10 show unsatisfactory fits, while Figure 9-11 shows a satisfactory fit for the decays "A" and "B".

20-11 When you are happy with the fit you can print the results directly from within the dialogue box (Print), or you can send the data to the clipboard, either as text or as image (**Copy as Text**, **Copy as Image**), for use in other software applications. To copy the fit plot, close the tail fit dialogue box, then press **Ctrl+C** on the graph.

Fit results (file suffix **F1**) and residuals (file suffix **F1R**) may be reviewed any time later by looking up the properties of the files with suffices **F1** and **F1R**, respectively. The properties of the residuals file also contain extended fit parameters, such as the Durbin-Watson parameter (Refer to chapter 9.4.4 page 109 for details).

During fitting, two potential error messages may be displayed:

- (1) "A matrix singularity occurred during fitting": One of your estimated lifetime guesses was too small, enter a value that is bigger and repeat the fit.
- (2) "A suitable fit was not found after 500 iterations": The fitting process was not convergent and was stopped after 500 iterations. This error occurs typically when trying to fit complex decays and the initial guesses are not sufficiently close to the "true" values. Another problem may be a mismatch between the decay and the IRF, i.e. the measured IRF does not represent the real function that was responsible for generating the measured sample decay.




Figure 9-7: Single Exponential Reconvolution Fit of example Decay "A"




Figure 9-8: Double Exponential Reconvolution Fit of example Decay "A"




Figure 9-9: Single Exponential Reconvolution Fit of sample Decay "B"




Figure 9-10: 2-exponential Reconvolution Fit of sample Decay "B"




Figure 9-11: 3-exponential Reconvolution Fit of sample Decay "B"

9.4. Fit Quality Parameters

The quality of the fit result can be evaluated in several ways. If the fit result is entirely wrong, then a simple visual comparison between the raw data and the fitted curve might be sufficient to find the reason. In most cases, however, the visual comparison between these two data is not sufficient. Other parameters need to be calculated to allow a much more precise fit evaluation.

The most common parameters are the following:

9.4.1. The Reduced Chi-Square

The numerical procedure behind the search for the best B_i and τ_i is the Marquardt-Levenberg algorithm. This is an iteration procedure which searches for the best B_i and τ_i by a controlled and directed minimisation of the "goodness of fit", χ_g^2 , which is defined as:

$$\chi_g^2 = \sum_k w_k^2 (F_k - S_k)^2$$

(k is the index for the individual data points to be fitted, the sum expands over all these data points.) The w_k are the weighting factors for the individual data points. Using the correct weighting factors for a specific set of raw data is important. The correct type of weighting factors is determined by the type of noise specific to the data and hence is inherited from the method which was used to collect the data. For example, lifetime data acquired by TCSPC or gated single photon counting (MCS) obey Poissonian noise statistics with the well defined weighting factor for each data point (S_k) of $w_k = 1/\sqrt{F_k}$. Data acquired by an oscilloscope obey Gaussian noise statistics with the $w_k = \text{const}$.

Dividing χ_g^2 by the number of free parameters n (which is approximately the number of fitted data points subtracted by the number of lifetime parameters used in the fit) will result in

$$\chi^2 = \sum_k w_k^2 \frac{[F_k - S_k]^2}{n}$$

χ^2 is called the "reduced chi-square"; it is the scaled "goodness of fit". The reduced chi-square has a distinct advantage over the goodness of fit in that its value is independent of the number of data points and the number of fitting parameters. This allows one to compare different fits.

For Poisson distributed data ($w_k = 1/\sqrt{S_k}$) the reduced chi-square has the theoretical limit 1.0. χ^2 -values above Unity indicate a bad fit result, although values of about 1.1, 1.2 or even 1.3 are acceptable under certain conditions. If the fitting range has been inappropriately chosen, χ^2 can be slightly less than 1.0.

Principally, one always needs to distinguish between the chi-square (χ_g^2) and the reduced chi-square (χ^2). As for fit evaluation only the reduced chi-square is used, the word "reduced" is omitted in this manual.

9.4.2. The Residual Data

Using the fit result data (F_k), the measured raw lifetime data (S_k), and the appropriate weighting factors (w_k - see next paragraph) the residual data can be calculated as follows:

$$X_k = w_k (F_k - S_k)$$

The residual data are the difference between the fitted curve and the original data to be fitted, weighted by the standard deviation of each data point. A good fit should give a residual curve that only contains random noise distributed around Zero. Any deviation from the randomness would give an indication for

a misfit, either because the appropriate exponential model contains more exponential terms than were used in the fit or because of instrumental artefacts.

9.4.3. The Autocorrelation Data

The autocorrelation function of the residuals is defined as

$$Z_k = \frac{\frac{1}{n_H - k - n_L} \sum_{i=n_L}^{n_H-k} X_i X_{i+k}}{\frac{1}{n_H - n_L} \sum_{i=n_L}^{n_H} X_i^2},$$

with n_L and n_H being the lower and upper limit of the fitting range, respectively. The residual autocorrelation data (Z_k) show more clearly than the residual data themselves (X_k) whether the residuals are fully randomly distributed or whether there is a repetitive pattern in the residuals. As each residual data point is correlated with itself the first autocorrelation data point is always 1.0. All other autocorrelation data points should be randomly distributed around Zero if the residuals are "clean".

The facility to calculate the autocorrelation function of the residuals will be available, if either a fit plot or a 2D display of a residual curve is in focus. If this is given, the autocorrelation function may be calculated via **Data > Autocorrelation....**

Upon opening of the dialogue box, simply click **Apply** to calculate the autocorrelation function. The function will be displayed in a new window and the file will carry the extension suffix **1RA**, **2RA**, etc, for single, double, etc fits.

9.4.4. The Durbin-Watson Parameter

This parameter is used to evaluate the quality of the fit. It is defined as

$$DW = \frac{\sum_{i=n_L+1}^{n_H} (X_i - X_{i-1})^2}{\sum_{i=n_L}^{n_H} X_i^2},$$

with n_L and n_H being the lower and upper limit of the fitting range, respectively. The Durbin Watson parameter can only be considered in absolute terms for a defined number of exponentials in the tried model. DW-values of less than 1.7, 1.75 and 1.8 are indicative for poor fits to single, double, and triple exponential decay models.

The Durbin-Watson Parameter is given in the table of properties of Time Scan Residuals (Figure 8-12)




Figure 9-12: Table of Properties of Time Scan Residuals

9.5. TRES Data Slicing

TRES-Map measurements (refer to section 6.5 on page 83) can be converted into time-resolved emission (or excitation) spectra by "slicing". The slices represent spectra of the sample's emission for given time windows. These spectra can provide valuable information as they show the dynamics of the emission during the sample's decay process.

Procedure 21: To convert a TRES-Map measurement into TRES...

- 21-1 Click **Data > TRES data slicing...**. This will open up a dialogue box as shown in Error! Reference source not found..
- 21-2 Define **Start time**, **Stop time**, and **No. of slices**. Update the **Current Slice Setting Information** by using the **Apply** button. This will not yet generate the data slices. You can control the width of the data slice by changing one of the three input parameters.
- 21-3 You may wish to set up each slice individually. You can do this by accessing the **advanced** setup.
- 21-4 Press **OK** to process the data slicing. A new window will be created, and the spectra within the new graph will be labelled with the starting time of each data slice.




Figure 9-13: TRES-Map with dialogue for Data Slicing




Figure 9-14: Time-Resolved Emission Spectra (TRES)

Formulas and Definitions

Definitions of Chromaticity	
Tristimulus functions	$X = \int_0^{\infty} I(\lambda) \bar{x}(\lambda) d\lambda$ $Y = \int_0^{\infty} I(\lambda) \bar{y}(\lambda) d\lambda$ $Z = \int_0^{\infty} I(\lambda) \bar{z}(\lambda) d\lambda$ (1-2)
CIE1931 Chromaticity functions x & y	$x = \frac{X}{X+Y+Z}$ $y = \frac{Y}{X+Y+Z}$ (3-4)
CIE 1976 Chromaticity functions u' and v'	$u' = \frac{4X}{X+15Y+3Z}$ $v' = \frac{9Y}{X+15Y+3Z}$ (5-6)
Exponential Decay Kinetics (example: 3-exponential model function)	
Exponential model function	$R(t) = B_1 \exp\left\{-\frac{t}{\tau_1}\right\} + B_2 \exp\left\{-\frac{t}{\tau_2}\right\} + B_3 \exp\left\{-\frac{t}{\tau_3}\right\}$ (7)
Relative concentration of the second component	$C_2 = \frac{B_2}{B_1+B_2+B_3}$ (8)
Relative fluorescence intensity of the second component, as a percentage	$\varphi_2 = \frac{B_2\tau_2}{B_1\tau_2+B_2\tau_2+B_3\tau_3} \cdot 100\%$ (9)
Intensity average lifetime of the entire fluorescence decay process	$\langle\tau\rangle_{int} = \frac{B_1\tau_1^2+B_2\tau_2^2+B_3\tau_3^2}{B_1\tau_1+B_2\tau_2+B_3\tau_3}$ (10)
Amplitude average lifetime of the entire fluorescence decay process	$\langle\tau\rangle_{amp} = \frac{B_1\tau_1}{B_1+B_2+B_3}$ (11)
Reconvolution of the Fluorescence Decay Model	
Exponential model function	$S(t) = \int_0^t E(t')R(t-t')dt'$ (12) <p><i>S(t)</i> – measured fluorescence decay <i>E(t)</i> – measured instrumental response function <i>R(t)</i> – theoretical sample decay model function</p>

Glossary

Acronym	Definition
CCD	Charge-Coupled Device
EHT	extra-high tension
EPL	Picosecond pulsed diode laser
EPLED	Picosecond pulsed LED
IRF	Instrument response function
MCP-PMT	Multi-Channel Plate Photomultiplier
MCS	Multi-Channel Scaling
NIR	Near Infrared
NIR-PMT	Near Infrared Photomultiplier
PMT	Photomultiplier
TCSPC	Time-Correlated Single Photon Counting
TRES	Time-Resolved Excitation/Emission Spectra
TRES-Map	Map of Time-Resolved Measurements for the purpose of generating TRES
UV	Ultraviolet
VUV	Vacuum Ultraviolet (radiation absorbed or scattered by air)

Index

Absorption	11, 46, 63
Analogue detector.....	28
Analysis Options	46
Anisotropy.....	11, 60
Attenuator	25
Bandwidth	33
CB1	4, 6, 7, 8, 9, 10
CCD.....	21
Combine.....	44
Correction File	46
Correction Scan	59
Crop Range	44
Cryostat	11
Cuvette.....	11, 17, 27
Czerny-Turner.....	18
Dark Count Rate	28, 29
Data	
Combine.....	44
Crop Range.....	44
Exp. Reconvolution Fit	45
Exp. Tail Fit.....	45
Normalise	44
Scale	31, 44
Subtract Baseline.....	44
Data Acquisition Card	4, 5, 6, 7, 8, 9, 10
Detector.....	3, 4, 5, 6, 7, 8, 9, 10, 11, 17, 18, 20, 21, 24, 25, 28, 29, 33
Detector Gating.....	39, 40, 42, 53, 68, 72, 109
Differentiate	44
Diode.....	13, 17
Dwell Time	15, 56
Em Detector	33
Em Wavelength.....	33
Emission Map.....	60
Emission Scan	53, 56, 70
Entrance Port	18, 19
EPL.....	13, 17
EPR Dewar.....	11
Excitation Scan.....	57
Excitation Wavelength.....	33
Exp. Reconvolution Fit.....	45
Exp. Tail Fit.....	45
FAST	97
File Icons.....	31
Filter.....	17, 24, 25
FLS1000.....	4, 5, 6, 7, 8, 9, 10, 13
Fluoracle	3, 30, 32, 34, 44, 45, 46, 67, 88, 97
Fluorescence	1, 3, 11, 18, 53, 109
Focal Length.....	18
Grating.....	14, 18, 19, 20, 21, 22, 26
InGaAs Detector	29
Instrument Response Function.....	45
Integrate	44
Integrating Sphere.....	27, 46
Intensity attenuator.....	33
Kinetic Scan	67, 70
Lamp Frequency	39, 40, 42, 53, 68, 72
Lanthanide.....	7
Laser.....	17
Lens Assembly.....	15, 26, 27
Lens Optics.....	26
Lifetime Measurement.....	1, 7, 14, 15, 30, 44, 45, 46, 67, 82
Kinetic.....	67, 70
Manual.....	72, 76
Multiple.....	80
Light Source	3, 4, 5, 6, 7, 8, 9, 10, 13, 17, 29, 33
Liquid Nitrogen.....	11, 29
Logarithmic Scale.....	31
Magnetic Stirrer.....	11, 27
Manual Lifetime Measurement	72, 76
Measurement Container.....	31
Measurement Container Icons.....	31
Measurement Options.....	46
Measurement Plot Icons.....	31
Monochromator	3, 14, 15, 18, 20, 21, 22, 24, 26
Double Monochromator.....	11
Emission Monochromator.....	11, 18, 20, 21
Excitation Monochromator	14, 15, 18, 21, 26
Multi-channel Scaling.....	3, 15
Multiple Lifetime Measurement.....	80
Multiple Scans.....	31, 32
Nanosecond	3, 5, 8, 9, 10, 19
Nanosecond Flashlamp	5, 8, 9, 10, 19
Near Infrared Detector	11
nF Lamp Setup	76

Nitrogen.....	11	Scale Icons.....	31
Normalise.....	42, 44	Scan Display Area	32
Options		Scan Icons.....	30
Analysis Options.....	46	Scan View Icons.....	31
Measurement Options.....	46	Setup	
Plot Defaults	31, 46	nF Lamp Setup.....	76
Sample Holder Options.....	27, 46	Signal Rate	30, 32, 33, 63
Ozone.....	14	µF Lamp Setup.....	72
Peak Value.....	90	Shutter	22
Phosphorescence	1, 3, 7, 11	Signal Rate	30, 32, 33, 54, 63
Photomultiplier	4, 5, 6, 7, 8, 9, 10, 28, 29	Silica.....	26
Gated	11	Single Photon Counting.....	3, 15, 53
Microchannel Plate.....	29	Smooth	44
Red Sensitive	28	Source.....	33
Photon.....	3, 5, 13, 18, 28, 56, 67	Spectral Measurement	30, 63
Picosecond.....	13, 17	Correction Scan	59
Picosecond Pulsed Diode Laser	13, 17	Emission Map.....	60
Plate Reader.....	11, 27	Emission Scan	56, 70
Plot Defaults.....	46	Excitation Scan.....	57
Plot options.....	42	Synchronous Map.....	62
Polarisation	22	Status Bar	32
Polariser.....	11	Steady State	7, 13, 28
Power Supply	29	Stray Light	11, 22, 24
Propagation Time	17	Subtract Baseline	44
Properties.....	31, 32, 57	Swing Mirror.....	18, 20, 21
Pulse Repetition Rate	15	Switch Lamp On/Off.....	76
Pulsed Light Emitting Diode.....	13, 17	Synchronous Map	62
Quantum Yield.....	28	Tail Fit.....	98, 101
Reconvolution Fit.....	101	TCC1.....	5, 8, 9, 10
Reference Detector	25	Temperature	11, 27, 29, 46
Repetition Rate	13, 15	Time Correlated Single Photon Counting (TCSPC)	3, 5, 13, 15, 29, 53
Sample Chamber.....	3, 17, 18, 20, 21, 24, 25, 26, 27, 63	Time Range.....	67
Sample Cooling.....	11, 27	Time Resolved.....	3, 7, 28, 29, 67
Sample Holder.....	11, 17, 27, 33, 46	Tool Bar	30
Sample Holder Options	27, 46	Ultraviolet.....	14, 16, 22, 25, 53
Sample Module	22, 30	Xenon Flash Lamp	15, 16
Sample Position.....	25, 26	Xenon Lamp.....	3, 4, 7, 9, 10, 13, 14, 16, 19
Sample Turret.....	27	Y-axis.....	31, 56
Samples.....	3, 11, 13, 15, 17, 18, 20, 21, 24, 25, 26, 27, 33, 45, 46, 63, 98	Zoom	31
Scale	31, 44	µF Lamp Setup.....	72

Table of Figures

Figure 1-1: General layout of, a) the FLS1000-SS-s and b) the FLS1000-DD-s.....	4
Figure 1-2: General layout of a) the FLS1000-xS-t with EPL/EPLED, b) the FLS1000-SS-t with nF920.....	5
Figure 1-3: General layout of an FLS1000-SS-m	6
Figure 1-4: General layout of an FLS1000-SS-sm.....	7
Figure 1-5: General layout of an FLS1000-SS-tm.....	8
Figure 1-6: General layout of an FLS1000-SS-st.....	9
Figure 1-7: General layout of a) an FLS1000-SS-stm, b) an FLS1000-DD-stm.....	10
Figure 2-1: Xe2 Xenon Lamp.....	13
Figure 2-2: Spectral Output of the 450 W Xenon Arc Lamp	14
Figure 2-3: μ F2 Xenon Flashlamp.....	15
Figure 2-4: Spectral Output of the Xenon Flashlamp.....	16
Figure 2-5: Temporal Output of μ F2 Flashlamp at 300 nm (blue), 450 nm (green) and 600 nm (red). ...	16
Figure 2-6: Edinburgh Instruments Picosecond Pulsed Diode Laser	17
Figure 2-7: TGM325-X Single-Grating Excitation Monochromator	18
Figure 2-8: Schematic of a TGM325-X Single-Grating Excitation Monochromator	19
Figure 2-9: DTGM325-X Double-Grating Excitation Monochromator	19
Figure 2-10: The TGM325-X Swing Mirror Assembly.....	20
Figure 2-11: TGM325-M Single-Grating Emission Monochromator.....	20
Figure 2-12: Schematic of Single-Grating Emission Monochromator	21
Figure 2-13: DTGM325-M Double Grating Emission Monochromator	21
Figure 2-14: Sample Chamber, shown with Optional Polariser and 2 nd Attenuator at the Laser Port..	24
Figure 2-15: Single Cuvette Sample Holder with Lens Assembly, Light Cones and Filters.....	27
<i>Figure 2-16: Detector Assembly, comprising Cooled Head to be connected to the CO1 Power Supply</i>	28
Figure 2-17: Quantum Efficiency of Red-Sensitive Photomultiplier	29
Figure 3-1: Signal Rate Window – principal layout.....	32
Figure 3-2: Excitation Scan, typical Display and Scan Properties	34
Figure 3-3: Emission Scan, typical Display and Scan Properties.....	35
Figure 3-4: Synchronous Scan, typical Display and Scan Properties	35
Figure 3-5: Example of an Emission Correction Scan, typical Display and Scan Properties.....	36
Figure 3-6: Time Scan, typical Display and Scan Properties	36
Figure 3-7: Fitted Time Scans (red) and Residual (green) with the respective detailed Properties.....	37

Figure 3-8: Setup of default 2D Plot Options	38
Figure 3-9: Example for 2D Display of Spectral data	38
Figure 3-10: Example for 2D Display of Time-Resolved data (IRF blue, measured decay red)	39
Figure 3-11: Example of a 3D display.....	40
Figure 3-12: Colour Map of Spectral data with three Cross-Hairs and Intensity Colour Scale.....	41
Figure 3-13: Contour Plot of Spectral data with three Cross-Hairs and Coloured Contour Lines	41
Figure 3-14: Contour Plot with filled Colour Map and overlaid (white) Contour Lines	41
Figure 3-15: Colour Map of time-resolved data with one cross-hair and intensity colour scale	42
Figure 3-16: Display of numerical values of a scan.....	43
Figure 4-1: Dialogue Box "Correction Files"	48
Figure 4-2: Dialogue Box "Window Options"	48
Figure 4-3: Dialogue Box "Measurement Options"	49
Figure 4-4: Options for Signal Rate display	50
Figure 4-5: Dialogue Box "Analysis Options"	50
Figure 4-6: Dialogue Box "Export to ASCII"	51
Figure 4-7: Dialogue Box "Sample Modules"	51
Figure 5-1: Signal Rate window for setting up spectral measurements.....	54
Figure 5-2: Signal Rate window with emission signal exceeding the recommended level	55
Figure 5-3: Emission Scan Setup dialogue box	56
Figure 5-4: Excitation Scan Setup dialogue box	58
Figure 5-5: Synchronous Scan setup dialogue box.....	59
Figure 5-6: Emission Map Setup.....	61
Figure 5-7: Emission Map Setup – Excitation	61
Figure 5-8: Emission Map Setup – Emission.....	62
Figure 5-9: Synchronous Map Setup	63
Figure 5-10: Multiple Scan Setup	64
Figure 6-1: Signal Rate window for setting up time-resolved measurements.....	71
Figure 6-2: Kinetic Scan Setup.....	72
Figure 6-3: μ F2 lamp setup dialogue box.....	74
Figure 6-4: Signal Rate dialogue of Time-Resolved Measurement by MCS.....	75
Figure 6-5: Dialogue box for a manual time-resolved measurement by MCS	76
Figure 6-6: EPL/EPLED setup dialog box for TCSPC	78
Figure 6-7: Signal Rate dialogue for TCSPC measurements	79

Figure 6-8: Dialogue box for a Manual Time-Resolved measurement by TCSPC.....	79
Figure 6-9: TCC setup: control over 5% rule for TCSPC measurements.....	80
Figure 6-10: Setup of Multiple Time-Resolved Measurements	82
Figure 6-11: Emission TRES-Map setup dialogue box	84
Figure 6-12. TCSPC measurements setup dialogues for Forward Mode.....	87
Figure 6-13. TCSPC measurements setup dialogues for Reverse Mode	87
Figure 7-1. Batch Measurements setup dialogue.....	89
Figure 7-2: Batch Measurement progress window.....	89
Figure 8-1: Add spectral scans.....	90
Figure 8-2: Append dialogue box	91
Figure 8-3: Dialogue box for scaling data.....	92
Figure 8-4: Normalise Data To Peak Value dialogue box	92
Figure 8-5: Normalise Data At Wavelength dialogue box.....	93
Figure 8-6: Subtract baseline	93
Figure 8-7: Crop Range dialogue box.....	94
Figure 8-8: Smooth dialogue box.....	94
Figure 8-9: Differentiate dialogue box.....	95
Figure 8-10: Integrate dialogue box.....	95
Figure 8-11: Assignment of Correction Files.....	96
Figure 8-12: Dialogue box for Manual Emission Correction.....	96
Figure 8-13: Analysis Options with option to enable CIE Chromaticity Calculation	97
<i>Figure 8-14: Chromaticity Plots displayed in CIE 1931 (top) and CIE 1976 (bottom)</i>	<i>98</i>
Figure 9-1: Tail fit regions, two different examples	101
Figure 9-2: Tail Fit dialogue box.....	102
Figure 9-3: satisfactory Tail Fit	103
Figure 9-4: unsatisfactory Tail Fit.....	103
Figure 9-5: Fit regions for Reconvolution Fit Analysis	104
Figure 9-6: Reconvolution Fit dialogue box.....	105
Figure 9-7: Single Exponential Reconvolution Fit of example Decay "A"	106
<i>Figure 9-8: Double Exponential Reconvolution Fit of example Decay "A"</i>	<i>106</i>
<i>Figure 9-9: Single Exponential Reconvolution Fit of sample Decay "B"</i>	<i>107</i>
Figure 9-10: 2-exponential Reconvolution Fit of sample Decay "B"	107
<i>Figure 9-11: 3-exponential Reconvolution Fit of sample Decay "B"</i>	<i>107</i>

Figure 9-12: Table of Properties of Time Scan Residuals.....110

Figure 9-13: TRES-Map with dialogue for Data Slicing.....110

Figure 9-14: Time-Resolved Emission Spectra (TRES).....111

Table of Procedures

<i>Procedure 1: Accessing the optical compartments inside the sample chamber</i>	25
<i>Procedure 2: Starting the spectrometer</i>	47
<i>Procedure 3: Shutting down the spectrometer</i>	52
<i>Procedure 4: To set-up and perform an emission scan</i>	56
<i>Procedure 5: To set-up and perform an excitation scan</i>	57
<i>Procedure 6: To set-up and perform a synchronous scan</i>	59
<i>Procedure 7: To set-up and perform an Emission Map measurement</i>	60
<i>Procedure 8: To set-up and perform a Synchronous Map measurement</i>	62
<i>Procedure 9: To set-up and perform a multiple spectral measurement:</i>	63
<i>Procedure 10: To set-up and perform a kinetic measurement</i>	72
<i>Procedure 11: To set-up and perform a time-resolved phosphorescence measurement</i>	75
<i>Procedure 12: To set-up and perform a measurement of the IRF</i>	77
<i>Procedure 13: To set-up and perform a time-resolved fluorescence measurement</i>	78
<i>Procedure 14: To set-up and perform a measurement of the IRF</i>	81
<i>Procedure 15: To set-up and perform a Multiple Time-Resolved Measurement</i>	82
<i>Procedure 16: To set-up and perform an Emission TRES Measurement:</i>	83
<i>Procedure 17: To set-up and perform Batch measurements</i>	88
<i>Procedure 18: To generate a chromaticity plot</i>	97
<i>Procedure 19: To perform a tail fit</i>	101
<i>Procedure 20: To perform a reconvolution fit</i>	104
<i>Procedure 21: To convert a TRES-Map measurement into TRES</i>	110

© Edinburgh Instruments Ltd
2 Bain Square, Kirkton Campus, Livingston, EH54 7DQ UK
Phone: +44 (0) 1506 425300 Fax: +44 (0) 1506 425320
Email: eai.service@edinst.com
<http://www.edinst.com>
